



Reach Global India Pvt. Ltd.

Technical Data Sheet (TDS)

TDS No.:

RGI/TDS/RSC

Rev. No.:

00

Effective Date:

Page No.

Page 1 of 1

Surgeon cap

PRODUCT CODE	RSC
RELATED STANDARDS	Surgeon Cap according to ISO13485 standards.
INTENDED USE	Surgeon Cap provide perfect fit with full head coverage and control cross contamination by hair, to doctors and nurses while conducting surgeries.
OPERATION CONDITIONS	NA
INSTRUCTION TO USE	Do not use if found damaged/contaminated Do not steam sterilize Dispose after use

PRODUCT SPECIFICATIONS:

Items	Specifications	
Bouffant Cap	Material used	Nonwoven fabric spun bond
	Color	Blue
	Thickness	0.15 mm
	GSM	25-30 gram
	Length of cap	12 inch
	weight	25 gram

Packing and Storage:

Packaging	100 Bouffant Caps per bag
Storage	Store in room temperature.
Expiry	Two years from the date of manufacturing.

Reach Global Pvt. Ltd.

Office No. 1&2, Wing A, Sai Raj Regency, Pimple Saudagar, Pune - 411 027, Maharashtra, India.

+91 20 2740 9487 | info@orohealthcare.com | Customer Care No: +91 93707 88850/52/55 | www.orohealthcare.com

CIN U 74900PN2012PTC145221 GST: 27AAGCR1308D1ZS