



HEALTHCARE
RANGE



Committed to build safe and better world for healthcare industry

Think sterilization and safety | Think ORO



ABOUT ORO



The veins of the ORO legacy trace back to 2009, when the distribution started from a small rented shop. This worthy lineage is well accepted across all the corners of India. Within a short span of time by 2013 ORO gained a popular name in the Dental Infection Control product range. ORO achieved accelerated growth and went on to become a preferred need in Dentistry with commitment, perseverance, and quality services. The beginning of 2015 marked yet another significant milestone, where ORO started its export to international markets. We started marking our presence in the International exhibition with a whopping response. We introduced Sure-endo to quantify customers demand with high-quality dental infection control

Carried the mantle forth by introducing modern, innovative ideas with smart services. Our products are now acknowledged globally in over 40+ countries. With a clear vision and focused thoughts, our team is ready to take the brand from the nascent stages to greater glory with our own manufacturing facility in India. Our enthusiastic and young team is all set to abreast modern, innovative ideas to ensure that our brand paces ahead adding value to our esteemed customers. With a highly competent team of professionals and their unified efforts, ORO believes to deliver sustainable success with a promise of persistent quality and availability.



Vision

While retaining our current reputation of "your infection guard" Our vision is to become 'leading manufacturer and supplier' of barrier packaging & and associated products for controlled area.

Mission

Our Sterilization packaging products, Protective clothing and supplementary products are targeted to cater the needs of 'Healthcare segment' for their CSSD, Production, Quality Control & Assurance sections. We aim to support their needs by offering suitable products manufactured at 'state-of-the art' facility with desire to maintain the 'best availability'

VALUES

Belief

Our strong belief is our people and we work by taking ownership to deliver our promises persistently



Unify

We encourage and challenge each other with respect to achieve our common objective with fun and by helping each other unified



Futuristic

We plan ourselves by targeting long term and sustainable success.



Understanding

We are open and have courage to try new and learn from our mistakes. We get happiness in developing partnership with our customers by rising understanding of small things to offer a solution





CONTENT

Sterilization Packaging Page

- | | |
|-------------------------------|------|
| 1. Medical Grade Reel & Pouch | 4-7 |
| 2. Tyvek Reel & Pouch | 8-10 |
| 3. SMMS Wrap Sheet | 11 |
| 4. Crepe Sheet | 12 |
| 5. Masking Tape | 13 |

Safety Product Page

- | | |
|----------------------------|-------|
| 1. Isolation Gown | 20 |
| 2. SMS Lab Coat & Jacket | 21-22 |
| 3. Face Mask | 23 |
| 4. Surgeon & Bouffant Cap | 24 |
| 5. Shoe Cover | 25 |
| 6. DuPont Tyvek® 400 | 26 |
| 7. DuPont Tyvek® 500 Xpert | 27 |
| 8. DuPont Tyvek® 600 Plus | 28 |

Sterility Monitoring Page

- | | |
|---------------------------------|-------|
| 1. Process Indicator Label | 14 |
| 2. Indicator Tape | 15 |
| 3. Bowie Dick Test Pack & Sheet | 16 |
| 4. Chemical Indicator | 17-19 |

STERILIZATION REEL AND POUCH

Speedy

Simple

Selection

Medical grade Sterilization reels & pouches allow sterilization of packaging contents while maintaining sterility until the package is opened. Ensuring delivery without any cross contamination


- Tear, fibre-free opening and aseptic presentation
- Triple layer sealing makes it easy to peel and with less force
- Wide range of standard sizes allows optimum choice
- Tinted transparent, reinforced multilayer co-polymer PET/ CPP film
- Excellent gas permeability and microbial barrier properties
- Suitable with Steam, EO sterilizers



Order Guide for ORO Sterilization Reel

Code	PRODUCT NAME	Packaging (Reels/ctn)
RRM50	ORO STERILIZATION REEL 50MM X 200M	10
RRM55	ORO STERILIZATION REEL 55MM X 200M	10
RRM75	ORO STERILIZATION REEL 75MM X 200M	8
RRM100	ORO STERILIZATION REEL 100MM X 200M	6
RRM125	ORO STERILIZATION REEL 125MM X 200M	8
RRM150	ORO STERILIZATION REEL 150MM X 200M	4
RRM200	ORO STERILIZATION REEL 200MM X 200M	4
RRM250	ORO STERILIZATION REEL 250MM X 200M	4
RRM300	ORO STERILIZATION REEL 300MM X 200M	2
RRM350	ORO STERILIZATION REEL 350MM X 200M	2
RRM400	ORO STERILIZATION REEL 400MM X 200M	2
RRM450	ORO STERILIZATION REEL 450MM X 200M	2
RRM500	ORO STERILIZATION REEL 500MM X 200M	2
RRM550	ORO STERILIZATION REEL 550MM X 200M	2
RRM600	ORO STERILIZATION REEL 600MM X 200M	2

*Customized sizes available on demand




Oredr Guide for ORO Sterilization Reel Gusseted

Code	PRODUCT NAME	Packaging (Reels/ctn)
RRM50	STERILIZATION REEL GUSSETED 50MM X 100M	10
RRM55	STERILIZATION REEL GUSSETED 55MM X 100M	10
RRM75	STERILIZATION REEL GUSSETED 75MM X 100M	8
RRM100	STERILIZATION REEL GUSSETED 100MM X 100M	6
RRM125	STERILIZATION REEL GUSSETED 125MM X 100M	8
RRM150	STERILIZATION REEL GUSSETED 150MM X 100M	4
RRM200	STERILIZATION REEL GUSSETED 200MM X 100M	4
RRM250	STERILIZATION REEL GUSSETED 250MM X 100M	4
RRM300	STERILIZATION REEL GUSSETED 300MM X 100M	2
RRM350	STERILIZATION REEL GUSSETED 350MM X 100M	2
RRM400	STERILIZATION REEL GUSSETED 400MM X 100M	2
RRM450	STERILIZATION REEL GUSSETED 450MM X 100M	2
RRM500	STERILIZATION REEL GUSSETED 500MM X 100M	2
RRM550	STERILIZATION REEL GUSSETED 550MM X 100M	2
RRM600	STERILIZATION REEL GUSSETED 600MM X 100M	2

***Customized sizes available on demand**



STERILIZATION REEL AND POUCH APPLICATION GUIDE



Order Guide For Self Seal Sterilization Pouch

Code	PRODUCT NAME	Packaging
RPMS50X228	ORO STERILIZATION POUCH SELF-SEAL 50 X 228MM	200 pcs/bag
RPMS57X100	ORO STERILIZATION POUCH SELF-SEAL 57 X 100 MM	200 pcs/bag
RPMS57X125	ORO STERILIZATION POUCH SELF-SEAL 57 X 125 MM	200 pcs/bag
RPMS60X130	ORO STERILIZATION POUCH SELF-SEAL 60 X 130MM	200 pcs/bag
RPMS60X363	ORO STERILIZATION POUCH SELF-SEAL 60 X 363MM	200 pcs/bag
RPMS70X255	ORO STERILIZATION POUCH SELF-SEAL 70 X 255MM	200 pcs/bag
RPMS75X280	ORO STERILIZATION POUCH SELF-SEAL 75 X 280MM	200 pcs/bag
RPMS90X165	ORO STERILIZATION POUCH SELF-SEAL 90 X 165MM	200 pcs/bag
RPMS90X260	ORO STERILIZATION POUCH SELF-SEAL 90 X 260MM	200 pcs/bag
RPMS100X125	ORO STERILIZATION POUCH SELF-SEAL 100 X 125MM	200 pcs/bag
RPMS100X150	ORO STERILIZATION POUCH SELF-SEAL 100 X 150MM	200 pcs/bag

***Customized sizes available on demand**



Order Guide For Heat Seal Sterilization Pouch

Code	PRODUCT NAME	Packaging
RPM50X203	ORO STERILIZATION POUCH HEAT-SEAL 50 X 203MM	200 pcs/bag
RPM57X70	ORO STERILIZATION POUCH HEAT-SEAL 57 X 70MM	200 pcs/bag
RPM57X95	ORO STERILIZATION POUCH HEAT-SEAL 57 X 95MM	200 pcs/bag
RPM60X103	ORO STERILIZATION POUCH HEAT-SEAL 60 X 103MM	200 pcs/bag
RPM60X337	ORO STERILIZATION POUCH HEAT-SEAL 60 X 337MM	200 pcs/bag
RPM70X230	ORO STERILIZATION POUCH HEAT-SEAL 70 X 230MM	200 pcs/bag
RPM75X254	ORO STERILIZATION POUCH HEAT-SEAL 75 X 254MM	200 pcs/bag
RPM90X138	ORO STERILIZATION POUCH HEAT-SEAL 90 X 138MM	200 pcs/bag
RPM90X233	ORO STERILIZATION POUCH HEAT-SEAL 90 X 233MM	200 pcs/bag
RPM100X103	ORO STERILIZATION POUCH HEAT-SEAL 100 X 103MM	200 pcs/bag

***Customized sizes available on demand**



TYVEK REEL & POUCH

Speedy

Simple

Selection

Tyvek is a nonwoven PE fabric which is a trademark product of Dupont. Tyvek is widely preferred for its outstanding strength, durability & tear resistant

- Excellent microbial barrier properties against germs, bacteria and moisture
- Extremely high resistance to penetration and fibre-free peeling
- Easy peel, seal and aseptic presentation
- Clear and accurate color change
- Long Shelf Life which can maintain sterility for at least five years if package integrity is not compromised



Order Guide for Sterilization Reel DuPont Tyvek® 2FS Plasma 70M


Code	PRODUCT NAME	Packaging
RRT2FSP5070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 50MM X 70M	10
RRT2FSP5570	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 55MM X 70M	10
RRT2FSP7570	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 75MM X 70M	8
RRT2FSP10070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 100MM X 70M	6
RRT2FSP12570	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 125MM X 70M	8
RRT2FSP15070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 150MM X 70M	4
RRT2FSP20070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 200MM X 70M	4
RRT2FSP25070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 250MM X 70M	4



Order Guide for Sterilization Reel DuPont Tyvek® 2FS Plasma 70M

Code	PRODUCT NAME	Packaging
RRT2FSP30070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 300MM X 70M	2
RRT2FSP35070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 350MM X 70M	2
RRT2FSP40070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 400MM X 70M	2
RRT2FSP45070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 450MM X 70M	2
RRT2FSP50070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 500MM X 70M	2
RRT2FSP55070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 550MM X 70M	2
RRT2FSP60070	ORO STERILIZATION REEL DUPONT TYVEK 2FS PLASMA 600MM X 70M	2

STERILIZATION REEL AND POUCH APPLICATION GUIDE





TYVEK REEL & POUCH

Order Guide for Sterilization Pouch Self-Seal DuPont Tyvek® 2FS Plasma

Code	Product Name	Packaging
RPTS2FSP60X288	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 60 X 228MM	200 Pcs / Box
RPTS2FSP57X100	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 57 X 100MM	200 Pcs / Box
RPTS2FSP57X125	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 57 X 125MM	200 Pcs / Box
RPTS2FSP60X130	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 60 X 130MM	200 Pcs / Box
RPTS2FSP60X363	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 60 X 363MM	200 Pcs / Box
RPTS2FSP70X255	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 70 X 255MM	200 Pcs / Box
RPTS2FSP75X280	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 75 X 280MM	200 Pcs / Box
RPTS2FSP90X165	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 90 X 165MM	200 Pcs / Box
RPTS2FSP90X260	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 90 X 260MM	200 Pcs / Box
RPTS2FSP100X125	STERILIZATION POUCH SELF-SEAL DUPONT TYVEK 2FS PLASMA 100 X 125MM	200 Pcs / Box

Order Guide for Sterilization Pouch Heat-Seal DuPont Tyvek® 2FS Plasma

Code	Product Name	Packaging
RPTS2FSP60X203	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 60 X 203MM	200 Pcs / Box
RPTS2FSP57X70	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 57 X 70MM	200 Pcs / Box
RPTS2FSP57X95	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 57 X 95MM	200 Pcs / Box
RPTS2FSP60X103	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 60 X 103MM	200 Pcs / Box
RPTS2FSP60X337	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 60 X 337MM	200 Pcs / Box
RPTS2FSP70X230	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 70 X 230MM	200 Pcs / Box
RPTS2FSP75X254	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 75 X 254MM	200 Pcs / Box
RPTS2FSP90X138	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 90 X 138MM	200 Pcs / Box
RPTS2FS90X233	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 90 X 233MM	200 Pcs / Box
RPTS2FSP100X103	STERILIZATION POUCH HEAT-SEAL DUPONT TYVEK 2FS PLASMA 100 X 103MM	200 Pcs / Box




SMS & SMMS WRAP SHEET

Simple

Convenient

Affordable

SMS & SMMS are made with layers of spunbound & meltblown overlapping each other



- Single use sheets for suitability for various standard & user define usage
- Sturdy material provides convenience in handling bulk product
- Excellent water-repellent material
- Smart Packaging for convenient usage
- Offers high tensile strength



*colour available on demand

GSM	CODE	PACKING
35 GSM	RSWSS3550x50B	700 pcs/ctn
	RSWSS3560x60B	700 pcs/ctn
	RSWSS3575x75B	400 pcs/ctn
	RSWSS3580x80B	400 pcs/ctn
	RSWSS3590x90B	400 pcs/ctn
	RSWSS35100x100B	400 pcs/ctn
	RSWSS35120x120B	300 pcs/ctn

GSM	CODE	PACKING
43 GSM	RSWSS4350x50B	600 pcs/ctn
	RSWSS4360x60B	600 pcs/ctn
	RSWSS4375x75B	350 pcs/ctn
	RSWSS4380x80B	350 pcs/ctn
	RSWSS4390x90B	400 pcs/ctn
	RSWSS43100x100B	400 pcs/ctn
	RSWSS43120x120B	250 pcs/ctn

GSM	CODE	PACKING
60 GSM	RSWSS6050x50B	500 pcs/ctn
	RSWSS6060x60B	500 pcs/ctn
	RSWSS6075x75B	250 pcs/ctn
	RSWSS6080x80B	300 pcs/ctn
	RSWSS6090x90B	300 pcs/ctn
	RSWSS60100x100B	300 pcs/ctn
	RSWSS60120x120B	200 pcs/ctn

*Customized sizes available on demand



CREPE SHEET

Reliable

Secure

Sustainable

Crepe sheets are processed from the direct extract of tree pulp and are made of 100% natural cellulose fibers, providing an excellent porosity & barrier for bacterial permeability.

- Ecofriendly trusted packaging material
- Soft, particle free fabrics provide most secure packaging
- Excellent Porosity providing high bacterial barrier
- High-Quality light weight material
- Smart packaging for convenient usage



*Colour available

GSM	CODE	PACKING
45 GSM	RSWSC4540x40G	1300 Pcs/ctn
	RSWSC4550x50G	1300 Pcs/ctn
	RSWSC4560x60G	1200 Pcs/ctn
	RSWSC4575x75G	600 Pcs/ctn
	RSWSC4580x80G	600 Pcs/ctn
	RSWSC4590x90G	250 Pcs/ctn
	RSWSC45100x100G	250 Pcs/ctn
	RSWSC45120x120G	250 Pcs/ctn
	RSWSC45150x150G	200 Pcs/ctn

GSM	CODE	CODE
60 GSM	RSWSC6040x40G	1000 Pcs/ctn
	RSWSC6050x50G	1000 Pcs/ctn
	RSWSC6060x60G	1000 Pcs/ctn
	RSWSC6075x75G	500 Pcs/ctn
	RSWSC6080x80G	250 Pcs/ctn
	RSWSC6090x90G	250 Pcs/ctn
	RSWSC60100x100G	250 Pcs/ctn
	RSWSC60120x120G	200 Pcs/ctn
	RSWSC60150x150G	200 Pcs/ctn

*Customized sizes available on demand



MASKING TAPE

Strength

Durability

Comfort

These high-quality adhesive tapes are intended to be used outside the article pack or container intended to be sterilized for keeping the package integrity intact.

Made of saturated crepe paper with pressure sensitive adhesive provide high strength.

Excellent durability with high working temperature upto 134°C.

Soft comfortable tearing makes it easy to handle by the user during bulky packaging



CODE	PACKAGING	SIZE
RMT	48 pcs/case	19mm*50m



PROCESS INDICATOR LABEL

Ensure Precise Sterilization Tracking With Us

Process Indicator Labels are Type 1 Indicators which on Exposure to the sterilization process; Help differentiate between processed and unprocessed loads

- 100% Lead And Heavy Metal Free Indicator Inks
- Precise & Sharp Process Indicator Colour Change
- High Quality Non Toxic Adhesive
- Packaging 1 roll x 500
- ISO 11140-1 compliant



CODE	TYPE	PACKAGING
RPILS	Steam	500 labels / roll
RPILE	EO	500 labels / roll
RPILP	Plasma	500 labels / roll



INDICATOR TAPE

For Simple packaging & Easy Identification

The indicator tapes are handy to identify processed & unprocessed pack at a glance

- High quality non-toxic adhesive
- Compatible to use over paper & film both
- Precise and sharp process indicator color change
- Easy recognizable indicator stripes
- 100% lead and heavy metal free indicator inks
- CONFORMS ISO 11140-1 & EN867



BEFORE

AFTER

STEAM INDICATOR TAPE

CODE: RSMSIT
SIZE: 19MM x 50M
PACK SIZE: 48PCS/CTN



EO INDICATOR TAPE

CODE: RSMEIT
SIZE: 19MM x 50M
PACK SIZE: 48PCS/CTN



PLASMA INDICATOR TAPE

CODE: RSMFIT
SIZE: 19MM x 50M
PACK SIZE: 48PCS/CTN





BOWIE DICK TEST PACK & SHEET

Everyday monitoring... make it simple

STERILITY MONITORING

Everyday monitoring Bowie & Dick Test Pack and sheets are best suited for equipment malfunction check

- Detects autoclave equipment malfunctions
- Precise and sharp process indicator colour change provides distinctive pass/fail result
- Smart & compact packaging
- 100% lead and heavy metal free indicator inks
- **Self Adhesive Test Sheet**



ORO BOWIE DICK TEST PACK TYPE-2 121°C-15MIN /134°C-3.5MIN

CODE: RSMBD-121/15
 CONFORMS ISO11140-4:2007
 CYCLE TIME: 121°C-15MIN/134°C-3.5MIN
 PACK SIZE: 30PKT/BOX

BEFORE



Pass



AFTER

ORO BOWIE DICK TEST PACK TYPE-2 121°C-10MIN /134°C-3.5MIN

CODE: RSMBD-121/10
 CONFORMS ISO11140-4
 CYCLE TIME: 121°C-10MIN/134°C-3.5MIN
 PACK SIZE: 30PKT/BOX



Pass



ORO BOWIE DICK TEST PACK TYPE-2 121°C-8.3MIN /134°C-3.5MIN

CODE: RSMBD-121/8.3
 CONFORMS ISO11140-4
 CYCLE TIME: 121°C-11MIN/134°C-3.5MIN
 PACK SIZE: 30PKT/BOX



Pass



BOWIE DICK TEST SHEETS

CODE: RSMBDS
 CONFORMS ISO11140-4
 CYCLE TIME: 121°C-15MIN/134°C-3.5MIN
 PACK SIZE: 50 SHEETS/BAG



Pass





CHEMICAL INDICATOR

Precise and sharp... makes sterilization easy

Chemical indicators are the devices which helps in confirming the meet of all sterilization parameters confirming a flaw proof sterilization.

- Monitors critical cycle parameters of steam sterilization
- Precise and sharp process indicator color change
- Laminated for durability
- Smart & compact packaging
- 100% lead and heavy metal free indicator inks

CHEMICAL INDICATOR STRIP FOR STEAM

CHEMICAL INDICATOR TYPE-4

CODE: RSMCIT4-121/15
 CONFORMS ISO 11140-1
 CYCLE TIME: 121°C-15min/134°C-3.5min
 PACK SIZE-200 STRIPS/BOX


CODE: RSMCIT4-121/20
 CONFORMS ISO 11140-1
 CYCLE TIME: 121° C-20 min/134° C-7 min
 PACK SIZE-200 STRIPS/BOX



CHEMICAL INDICATOR TYPE-5

CODE: RSMCIT5-121/15
 CONFORMS ISO 11140-1
 CYCLE TIME: 121° C-15 min/134° C-3.5 min
 PACK SIZE-200 STRIPS/BOX

CODE: RSMCIT5-121/20
 CONFORMS ISO 11140-1
 CYCLE TIME: 121° C-20 min/134° C-7 min
 PACK SIZE-200 STRIPS/BOX





STERILITY MONITORING

CHEMICAL INDICATOR

Precise and sharp... makes sterilization easy

CHEMICAL INDICATOR TYPE-6

CODE: RSMCIT6-121/15

CONFORMS ISO 11140-1

CYCLE TIME: 121° C-15 min/134° C-3.5 min

PACK SIZE-200 STRIPS/BOX

CODE: RSMCIT6-121/20

CONFORMS ISO 11140-1

CYCLE TIME: 121° C-20 min/134° C-7 min

PACK SIZE-200 STRIPS/BOX



CHEMICAL INDICATOR STRIP FOR EO

CHEMICAL INDICATOR TYPE-4 EO

CODE: RSMCIE4

CONFORMS ISO 11140-1

PACK SIZE-200 STRIPS/BOX





STERILITY MONITORING

CHEMICAL INDICATOR

Precise and sharp... makes sterilization easy

CHEMICAL INDICATOR TYPE-5 EO

CODE: RSMCIE5

CONFORMS ISO 11140-1

PACK SIZE-200 STRIPS/BOX



CHEMICAL INDICATOR TYPE-6 EO

CODE: RSMCIE6

CONFORMS ISO 11140-1

PACK SIZE-200 STRIPS/BOX





SURGEON & ISOLATION GOWN

Protection

Safe

Comfort

The non-woven isolation gowns are highly breathable providing best protection

- Excellent fluid resistant nonwoven fabric
- High microbial barrier
- Great breathability giving all day comfort
- Ultimate protection against dirt & dry particulates
- Ideal for sterile area usage



Available GSM : 35,43 & 60

Code	Size	Type	packaging
RIGR	S, M, L, XL, XXL	REINFORCED STERILE	10 Pcs / bag
RIGE	S, M, L, XL, XXL	NON-STERILE -ELASTIC	10 Pcs / bag
RIGK	S, M, L, XL, XXL	NON-STERILE-KNITTED	10 Pcs / bag
RIGS	S, M, L, XL, XXL	STERILE SURGEON GOWN	10 Pcs / bag

***Customized size & GSM available on demand**



SMS LAB COAT & JACKET

Protection

Safe

Comfort

These are elegant professional wear for healthcare professionals

LAB COAT: Knee Length

JACKET: Hip Length

- Smart Look
- Knitted cuff & collar for comfort and complete enclosure
- Soft breathable & fluid resistant fabric assuring all day protection
- Pockets & snap closure





SMS LAB COAT

Code	Size	Packaging
RLCSB	S, M, L, XL, XXL	10 pcs/bag

SMS JACKET

Code	Size	Packaging
RJSB	S, M, L	10 pcs/bag
RJSL	S, M, L	10 pcs/bag
RJSP	S, M, L	10 pcs/bag

* Please check for MOQ for colors except blue



* colour available on demand



FACE MASK

Protection

Safe

Comfort

ORO non absorbant & non-irritating masks remains cool on the skin providing ultimate comfort

- Fluid Resistant mask
- Adjustable nose strip
- Soft earloop for comfortable wear
- TYPE IIR mask conforming to EN14683:2019+AC:2019
- $\geq 98.4\%$ BFE, Superior Meltblown
- Earloop & Tie-on variants available
- Available in 3PLY & 4PLY

*colour available for 3 PLY & 4 PLY



Code	Product Name	Packaging
RFM3E	Face Mask 3 Ply Ear Loop (Blue, Black, Pink, White)	50 Pcs / box
RFM4E	Face Mask 4 Ply Ear Loop (Blue, Black Pink, White)	50 Pcs / box
RFM3T	Face Mask 3 Ply Tie On (Blue, Black Pink, White)	50 Pcs / box
RFM4T	Face Mask 3 Ply Tie On (Blue, Black Pink, White)	50 Pcs / box



SURGEON & BOUFFANT CAP

Protection

Safe

Comfort

The disposable surgical caps are meant to be used by doctors and nurses while conducting surgeries. ORO caps are manufactured using optimum quality fabric, which gives great comfort & highly breathable. Soft fabric with smooth texture makes the caps light weighted.

Bouffant caps provide the perfect fit with full head coverage, and control cross-contamination by hair. These caps are high breathability, light weight and antistatic, with an elegant look.



Code	Product Name	Type	Packaging
RSCNW	SURGEON CAP	Blue	100 pcs/bag
RBCS18B, RBCS21B	BOUFFANT CAP BLUE	Strip 18" & 21"	100 pcs/bag
RBCS18W, RBCS21W	BOUFFANT CAP WHITE	Strip 18" & 21" White	100 pcs/bag



SHOE COVER

Protection

Safe

Comfort

The ORO universal size Shoe Cover made of Non-Woven or LDPE Fabric. It is water resistant, anti-static & Anti-Skid. The shoes gives great comfort for long hours working with great comfort & extreme hygiene.



Code	Product Name	Packaging
RSHCN	SHOE COVER NON-WOVEN	100 Pcs / bag
RSHCL	SHOE COVER LDPE	100 Pcs / bag



DuPont Tyvek® 400

FEATURES

- Hooded coverall
- Elastic opening for tighter fit at wrist and ankle
- Longer zipper extends to chin for complete coverage of neck area
- Storm flaps



Antibacterial



Breathable



Waterproof



Product Code	RGTYC400
Colour	White
Sizes	S, M, L, XL, XXL

Tyvek® is a registered trademark of DuPont™.



DuPont™ Tyvek® 500 Xpert

FEATURES

- Hooded coveralls with elastic
- Ergonomic- protective design
- Stitched external seams
- Elasticated waist, ankles and face
- Chemical protective clothing, Category III, Type 5-B and 6B
- Anti-static treatment (EN 1149-5) on both sides



Product Code	RGTYC500
Colour	White
Sizes	S, M, L, XL, XXL

Tyvek® is a registered trademark of DuPont™.



DuPont™ Tyvek® 600 Plus

FEATURES

- Hooded coveralls
- Stitched and over taped seams
- Thumb loops
- Self-adhesive zipper and chip flap
- Ergo-nomic - protective design
- Stitched external seams
- Elasticated waist, Ankles, face & wrist
- Chemical protective clothing, Category III, Type 4-B ,5-B and 6-B
- Anti-static treatment (EN 1149-5) on both sides



Product Code	RGTYC600
Colour	White
Sizes	S, M, L, XL, XXL

Tyvek® is a registered trademark of DuPont™.



ACCREDITATIONS





Reach Global India Pvt. Ltd.

- 📍 Reach Global India Private Limited.
Office No. 1&2, Survey No. 64/1, Wing-A, First Floor,
Sai Raj Regency, Pimple Saudagar, Pune: 411027
- 📍 Plot No I-22, Gat No. 1064, Chordia Industrial Park,
Naigaon, Shirwal, Taluka-Khandala,
Dist. Satara, State-Maharashtra (India)
Pin Code - 412801

📞 +91 20 2740 0987

📞 +91 93707 88850

✉️ info@orohealthcare.com

🌐 www.orohealthcare.com