



DENTAL CONSUMABLES & STERILIZATION



COMMITTED TO CREATING NO INFECTION ZONES

Today, the dental industry is growing at a rapid pace. With growing wealth and awareness of dental healthcare, customers are seeking treatment from clinics and outlets that employ high-end dental technologies, use best-in-class dental consumable products and maintain strict standards of hygiene.

ORO is an established and reputed player in dental consumable products with superior quality and innovative technology. A trusted brand with the most comprehensive range of products to promote exceptional level of hygiene in dentistry. ORO goes beyond the core values of safety and effectiveness. ORO's vision and initiative is to enable dental practices across India and abroad stay hygienic and completely devoid of any infections. ORO – not just in the business of dental consumables, but in the business of creating a No Infection Zone.





REACH GLOBAL

Redefining standards in the Healthcare Industry

Reach Global is a highly reputed name in the dental industry in terms of product quality and financial performance with best-in-class service. Over the years, Reach Global has consistently kept itself abreast of the latest technologies and stayed ahead of the curve in the industry.

Operating as a true one-stop solution, Reach Global has always provided its customer network a basket of high-range performance products and solutions at an optimum cost. Besides demonstrated performance, durability, low maintenance and safe operations, the brand has earned itself the niche of also being rated high on cost-efficiency. Thereby delivering optimum value with thorough professionalism in business operations.

With the principal focus on providing solutions which are patient, doctor & environment friendly, the brand's innovative products meet the highest standards in the dental industry.

Besides India, the representative office is also strategically located at Shenzhen, Southeast China, to bring swiftness and convenience in routine logistics and delivery operations.

With a commitment to deliver 'No Infection Zone', Reach Global is pillared on the values of high quality, transparent dealings and flawless distribution to stay at the forefront of meeting each and every requirement of the dental industry in India and overseas.

REACH GLOBAL

INFECTION BARRIERS

Barrier Film.....05
 Barrier Envelope.....05
 Sleeves.....06

SAFE STERILIZATIONS

Bowie dick Kit.....06
 Indicator Tapes.....07
 Safe Seal Pouches.....07
 Sterilization Rolls.....08
 Steam Indicator Strip.....08

IMPRESSION ACCESSORIES

Bite Tray.....09
 Fluoride Tray.....09
 Impression Tray.....10
 Mixing Bowl.....10
 Spatula.....11
 Dispensing Gun.....11
 Mixing Tips.....12

GERM-FREE EVACUATIONS

Intra Oral Tip.....13
 Saliva Ejector.....13
 Surgical Aspirator Tip.....14
 Vented Oral Evacuation Tip.....14

HYGIENIC DISPOSABLES

Dental Kit.....15
 Disposable Air / Water
 Syringe Tips.....15
 Micro Applicator.....16

SECURE DISPENSERS / ORGANIZERS

CompuLe Gun.....16
 Dental Plastic Tray.....17
 Endo Block.....17
 Endo Boxes.....18

RETRACTORS AND EXPANDERS

Cheek Retractor.....18
 Mouth Prob.....19

PROTECTIVES & ABSORPTIVES

Bouffant Cap.....19
 Face Mask.....20
 Face Shield.....20
 Implant Kit.....21
 Patient Bib.....21
 Dental Apron.....22

Protective Glass.....22

SMS Lab Coat & Jacket.....23
 Surgical Gown.....23
 Cotton Roll & Dispenser.....24
 Sponges.....24

SAFE KEEPERS

Bur holder.....25
 Denture Box.....25
 Retainer Box.....26

MISCELLANEOUS

Glass Dappen Dish.....26
 Prebent tips.....27
 Wooden Wedges.....27
 Prophy Cup & Brush.....28
 X-ray film holder.....28

INFECTION BARRIERS

BARRIER FILM



These perforated sheets with non-sticky edge are effective to reduce cross contamination with single-use and covers anything you might touch during procedure. Easy to protect the areas which are difficult to clean.

Smart packaging with self-dispense is best for direct usage with easy tear off .

SPECIFICATIONS

Code	Size	Packing
Blue, RBFB	4" x 6", 1200 sheets	1 wrap / box, 8 boxes / case
Clear, RBFC	4" x 6", 1200 sheets	1 wrap / box, 8 boxes / case

INFECTION BARRIERS

BARRIER ENVELOPES



Barrier Envelope protects the phosphor storage plates from mechanical damage and information loss caused by ambient light. The Envelope keeps the phosphor storage plate clean while reducing cross contamination.

Peel and stick end provide a consistent seal and tear off.

SPECIFICATIONS

Code	Size	Packing
RBE0	Size 0	100 pcs / box, 100 boxes / case
RBE1	Size 1	100 pcs / box, 100 boxes / case
RBE2	Size 2	100 pcs / box, 100 boxes / case

INFECTION BARRIERS

SLEEVES



Sleeves are important to maintain hygienic workflow and are vital to operator asepsis, which play crucial role in infection control by reducing clean up time and avoid cross-contamination.

We are specialized to provide sleeves for all varieties and sizes of your requirement.

SPECIFICATIONS

Type	Code	Size	Packing
Sensor Sleeve	RSLIC10.2X2	RSLIC10.2X2	500 pcs / box, 36 boxes / case
Scaler Sleeve	RSLLC4.3X2.95	RSLLC4.3X2.95	500 pcs / box, 36 boxes / case
Intra Oral Camera Sleeve	INSLV	RSLSC1.4X7	500 pcs / box, 36 boxes / case
Light Cure Sleeve	LCTC	RSLS8X1.6	200 pcs / box, 42 boxes / case

SAFE STERILIZATIONS

BOWIE DICK



The test pack consists of a series of air removal and steam penetration barriers. Single use test packs are designed for validation in sterilization process in the sterilization process.

- Colour change from Blue to Black will be clear and stable, and will not disappear easily
- Manufactured in the absence of lead or other heavy metals
- Design to spot any malfunction / lack of efficacy in the pre-vacuum sterilizer

SPECIFICATIONS

Code	Packing
RSMBD	120 kit / case (121°C-15min, 134°C-3.5min)

SAFE STERILIZATIONS

INDICATOR TAPES



Autoclavable adhesive tapes are used to conform whether a specific temperature has been attended. These are water-based and non-toxic chemical indicator tapes with stripes. Extra adhesive is provided to fix on all types of paper & plastic.

SPECIFICATIONS

Type	Code	Size	Packing
STEAM	ASST	19mm x 50m	48 Pcs / case
EO	AEST	19mm x 50m	48 Pcs / case
Plasma	APST	19mm x 50m	48 Pcs / case

SAFE STERILIZATIONS

SAFE SEAL POUCHES



These well-designed sterilization pouches come with a triple-layered seal, delivering three levels of safety for easy, effective and infection-free sterilization. These pouches conform to the three parameters of time, temperature & pressure steam / EO valid up to 6 months. Easy monitoring can be achieved through internal and external colour change indicators.

SPECIFICATIONS

Code	Size	Packing
RRM57x100	2.2**8.9"(57mm*100mm)	1000 pcs / box, 50 boxes / case
RRM57X125	2.2**4.9"(57mm*125mm)	6000 pcs / box, 30 boxes / case
RRM90X160	3.5**6.2"(90mm*160mm)	6000 pcs / box, 30 boxes / case
RRM70X255	2.7**10"(70mm*255mm)	4000 pcs / box, 20 boxes / case
RRM90X260	3.5**10"(90mm*260mm)	4000 pcs / box, 30 boxes / case
RRM100X280	3.9**11"(100mm*280mm)	4000 pcs / box, 20 boxes / case
RRM134X283	5.2**11"(134mm*283mm)	4000 pcs / box, 20 boxes / case
RRM190X360	7.4**14"(190mm*360mm)	2000 pcs / box, 10 boxes / case
RRM250X400	9.8**15" (250mm*400mm)	1000 pcs / box, 5 boxes / case
RRM300X400	11.8**15"(300mm*400mm)	1000 pcs / box, 5 boxes / case
RRM150x280	5.9**11.02 (150mm*280mm)	1000 pcs / box, 10 boxes / case
RRM330x450	12.9** 17.7" (330mm*450mm)	1000 pcs / box, 5 boxes / case

SAFE STERILIZATIONS

STERILIZATION ROLLS



Tripled sealed flat rolls are constructed from transparent PET/CPP multi-layers copolymer film and medical grade paper. The process indicators for steam and ethylene oxide are water-based, non-toxic inks which provide accurate colour change indication to confirm the process is achieved.

SPECIFICATIONS

Code	Size	Packing
RRM50	50mm(1")*200m	10 rolls / case
RRM55	55mm(2")*200m	10 rolls / case
RRM75	75mm(2.9")*200m	8 rolls / case
RRM100	100mm(3.9")*200m	6 rolls / case
RRM125	125mm(4.9")*200m	8 rolls / case
RRM150	150mm(5.9")*200m	4 rolls / case
RRM200	200mm(7.8")*200m	4 rolls / case
RRM250	250mm(9.8")*200m	4 rolls / case
RRM300	300mm(11")*200m	2 rolls / case
RRM350	350mm(13")*200m	2 rolls / case
RRM400	400mm(15")*200m	2 rolls / case

SAFE STERILIZATIONS

STEAM INDICATOR STRIP



Steam Indicator Strips are designed to be used in steam sterilizers operating at either 121°C for 15 minutes or 134°C for 3-5 minutes. The water-based non-toxic chemical indicator will show a clear and accurate colour change from blue to black after successful steam sterilization.

SPECIFICATIONS

Code	Type	Packaging
RSMCIT4	Class IV	200 pcs / box, 100 boxes / case
RSMCIT5	Class V	200 pcs / box, 100 boxes / case
RSMCIT6	Class VI	200 pcs / box, 100 boxes / case

IMPRESSION ACCESSORIES

BITE TRAY



Bite registration tray is an easy and economical way to reduce the risk of cross-contamination and infection. The robust construction of the single use dual-bite impression tray is made with polished plastic, smooth edges to ensure the patient's comfort. The net gives best holding capability to under-lock the material, and retains it in the same position.

SPECIFICATIONS

Code	Type	Packing
RBTAT	Anterior	35 pcs / bag, 20 bag / case
RBTAS	Assorted	50 pcs / bag, 20 bag / case
RBTFA	Full Arch	35 pcs / bag, 20 bag / case
RBTPT	Posterior	50 pcs / bag, 20 bag / case
RBTQD	Quadrant	30 pcs / bag, 20 bag / case
RBTSL	Sideless	35 pcs / bag, 20 bag / case

IMPRESSION ACCESSORIES

FLUORIDE TRAY




Disposable fluoride trays available in a variety of sizes are custom-fitted foam trays for fluoride gel, to control decay and strengthen teeth structure. They are designed to offer a natural, contoured fit and exceptional comfort for fluoride treatment.

SPECIFICATIONS

Code	Size & Colour	Packing
RFTL	Large, Blue	50 pcs / bag, 50 bags / case
RFTM	Medium, Yellow	50 pcs / bag, 60 bags / case
RFTS	Small, Green	50 pcs / bag, 80 bags / case

IMPRESSION ACCESSORIES

IMPRESSION TRAY



Perforated impression trays are made from strong, rigid materials to generate consistent impressions. These trays are distortion-free, pre-shaped, readily trimmed, and are flame-softened for adjustment. The handle features non-slip grip and are compatible with wide variety of impression materials. Available in disposable and autoclavable material up to 121°C.

SPECIFICATIONS

Code	Sizes	Packing
RIMTA1&2	Size 1&2	10 pcs / bag, 50 bags / case
RIMTA3&4	Size 3&4	10 pcs / bag, 50 bags / case
RIMTA5&6	Size 5&6	10 pcs / bag, 50 bags / case
RIMTA7&8	Size 7&8	10 pcs / bag, 50 bags / case
RIMTA9&10	Size 9&10	10 pcs / bag, 50 bags / case
RIMTA1-10	Size 1-10	10 pcs / bag, 50 bags / case

IMPRESSION ACCESSORIES

MIXING BOWL



Flexible PVC bowl work best for mixing impression material and is easy to re-use after cleaning. The material is suitable for both hard and soft mixtures.

SPECIFICATIONS

Code	Capacity	Colour	Packing
RMBL600	600ml	Blue	100 pcs / case
RMBM500	500ml	Green	100 pcs / case
RMBS280	280ml	Clear Green / Red	180 pcs / case
RMBXS150	150ml	Clear Green / Red	250 pcs / case



MIXING TIP



A wide range of mixing tip allows an easy choice. As per material suitability with dynamic proportion, and the collection suits for extra-light / light / regular body impression material, regular / heavy body impression material etc. The tips allow for easy, automatic mixing of impression material using the ORO dispensing gun, and are perfect to deliver convenient and precise material placement.

MIXING TIP

SPECIFICATIONS



Code	Type	Image	Application	Packing
RMTBFC1:1	Brown Flat Core 1:1		Dual Core build up material	50 pcs / bag, 50 bags / case
RMTBTRS1:1	Brown Tapered Resin Short 1:1		Resin Material	50 pcs / bag, 50 bags / case
RMTBTRL1:1	Brown Tapered Resin Long 1:1		Resin Material	50 pcs / bag, 50 bags / case
RMTYLBL1:1	Yellow Light Body Long 1:1		Light Body	50 pcs / bag, 50 bags / case
RMTYLSB1:1	Yellow Light Body Short 1:1		Light Body	50 pcs / bag, 50 bags / case
RMTGHBL1:1	Green Heavy Body Long 1:1		Heavy Body	50 pcs / bag, 50 bags / case
RMTGHBS1:1	Green Heavy Body Short 1:1		Heavy Body	50 pcs / bag, 50 bags / case
RMTPR/HB1:1	Pink Regular Regular/Heavy Body 1:1		Medium Body	50 pcs / bag, 50 bags / case
RMTTC&B1:10	Temporary Crown and Bridge 1:10		Crown and Bridge	50 pcs / bag, 50 bags / case
RMTBTC1:1	Blue Tapered Cement 1:1		Crown and Bridge	50 pcs / bag, 50 bags / case
RMTBTCL1:10	Brown Temporary Cement Long 1:10		Temporary Cement	50 pcs / bag, 50 bags / case
RMTBTCS1:10	Brown Temporary Cement Short 1:10		Temporary Cement	50 pcs / bag, 50 bags / case
RMTYLSB1:1	Blue Light Body Short 1:1		Light Body	50 pcs / bag, 50 bags / case
RMTRD	Red Dynamic		Red Dynamic	50 pcs / bag, 50 bags / case
RMTYD	Yellow Dynamic		380ml VPS	50 pcs / bag, 50 bags / case

IMPRESSION ACCESSORIES

SPATULA



A flat and smooth fore-end spatula is good for mixing. Autoclavable up to 121°C.

SPECIFICATIONS

Code	Colour	Packing
RSP	Pink, White, Yellow & Blue	50 pcs / bag, 10 bags / case

IMPRESSION ACCESSORIES

DISPENSING GUN



Dispensing gun works perfectly to dispense dental adhesive silicones and dental materials for perfect registration of dental impression. This sturdy device is made of high-quality plastic and consists of a release lever with a plunger, which coordinates simultaneously to load and remove the cartridges. Works perfectly with ORO mixing tips.

SPECIFICATIONS

Code	Type	Packing
RDG1:1X50ML	1:1/1:2 50 Mi	1 pc / box, 50 boxes / case
RDG1:1X25ML	1:1/1:4 25 Mi	1 pc / box, 50 boxes / case
RDG1:10X50ML	1:10/1:4 50 Mi	1 pc / box, 50 boxes / case

IMPRESSION ACCESSORIES

INTRA ORAL TIP



Oro Intra Oral Tips are dispensing tips ,especially designed for controlled release of core buildup material and other medium viscosity flowable materials.It ensures precise and accurate placement of materials even into the difficult to access areas.It is disposable ,therefore eliminates cross contamination from one patient to the other.The anti-slippery design prevents the tip from slipping or detaching from the automixing tips or cartridges to which they are attached.It is compatible with most existent mixing tips,therefore,ensures minimal wastage of material and easy clean-up after use.

SPECIFICATIONS

Code	Type	Packing
RIOTN1	N1	50 pcs/bag / 100 bags / case
RIOTN2	N2	50 pcs/bag / 100 bags /case
RIOTN3	N3	50 pcs/bag / 100 bags / case
RIOTN5	N5	50 pcs/bag / 100 bags / case
RIOTN8	N8	50 pcs/bag / 100 bags / case

GERM FREE EVACUATIONS

SALIVA EJECTOR



The non-sterile disposable saliva ejector is made of non-toxic material. The incorporated copper wire support to maintain a stable and precise shape. The non-removable and special shape tip guarantees optimal suction, preventing any damage to the Patient's oral mucosa.

SPECIFICATIONS

Code	Colour	Packing
RSEB	Blue	100 pcs / bag , 50 bags / case
RSEC	Clear	100 pcs / bag , 50 bags / case
RSEG	Green	100 pcs / bag , 50 bags / case

GERM FREE EVACUATIONS

SURGICAL ASPIRATOR TIP




ORO's surgical aspirator tips are made of high-quality polypropylene and designed to fit all high-volume evacuations. These tips are rigid enough to aspirate and retract the cheek without bending, and the 30° permanent angle allows better visual contact with the site, while the smooth edges enhance patient comfort. These can be disposed after use or autoclaved at 121°C.

SPECIFICATIONS

Code	Size & Colour	Packing
RSPT1/16B	1/16" Blue	25 pcs / bag, 50 bags / case
RSPT1/4G	1/4" Green	25 pcs / bag, 50 bags / case
RSPT1/8W	1/8" White	25 pcs / bag, 50 bags / case

GERM FREE EVACUATIONS

VENTED ORAL EVACUATION TIP



ORO's single use disposable tips contain smooth edges for maximum patient comfort. One end is vented and the other is non-vented. These tips fit all standard evacuator hose ends and are non-sterilized.

SPECIFICATIONS

Code	Size & Colour	Packing
RVOET6	6", white	100 pcs / bag, 10 bags / case
RVOET3	3", white	100 pcs / bag, 20 bags / case

HYGIENIC AND DISPOSABLES

DENTAL KIT



Dental examination kits are pre-sterile by EO and are ready-to-use. These kits are made with ABS material, and provide safety and convenience.

- Ready-to-use
- Minimises tray set-up time
- Hygienic (single-use)
- Eliminates cross-contamination issues

SPECIFICATIONS

Code	Size & Colour	Packing
DK - 001	Mirror, Probe & Tweezer Blue, White	3pcs / kit , 100 kits / box, 6 boxes / set

HYGIENIC AND DISPOSABLES

DISPOSABLE AIR / WATER SYRINGE TIPS



These tips have an inner core water passage channel that fits nicely into the syringe preventing water contamination in the air duct. The core design of the tips provides instant dry air and prevents cross-contamination while delivering superior air / water spray ratios.

SPECIFICATIONS

Code	Colour	Packing
RAWSDA	Assorted	250 pcs / box, 50 boxes / case
RAWSDB	Blue	50 pcs / box, 50 boxes / case
RAWSDC	Clear	250 pcs / box, 50 boxes / case

HYGIENIC AND DISPOSABLES

MICRO APPLICATOR



Micro applicator tip is non-linting and non-absorbent with a featured bendable handle. These tips are designed for precise placement of bonding agents, etchants, haemostatic solutions, cavity liners and disclosing solutions. They facilitate the proper placement of materials into difficult areas. Available in multiple colours for easy identification in a variety of sizes.

SPECIFICATIONS

Code	Colour	Packing
RMASFW	Superfine White	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMAFM	Fine Mauve	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMAFY	Fine Yellow	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMARB	Regular Blue	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMARG	Regular Green	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMARO	Regular Orange	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMARP	Regular Purple	100 pcs / tube, 4 tubes / box, 32 boxes / case
RMASBL	X - Slim, Black	100 pcs / tube, 4 tubes / box, 32 boxes / case

SECURE DISPENSER / ORGANIZER

COMPULE GUN



ORO's compule gun is designed to deliver the material with easy handling of the compule tips, i.e. direct and precise placement of the composite and compomer restorative into the cavity. These injection guns match most charging barrels, and are the best choice for sticky materials.

SPECIFICATIONS

Code	Size & Colour	Packing
CMG160	150mm x 100mm	160 pcs / case

SECURE DISPENSER / ORGANIZER

DENTAL PLASTIC TRAY




Experience the convenience of a plastic set up tray facility during any procedure. The trays are made from high quality plastic ensuring long durability. They have been meticulously designed to allow for systematic placement of instruments.

SPECIFICATIONS

Code	Type	Packing
RITA22x16.5, RITA33x24, RITA20x10	Autoclavable	160 pcs / case
RITD28x19B, RITD28x19G, RITD28x19O, RITD28x19W, RITD28x19Y	Disposable	160 pcs / case

SECURE DISPENSER / ORGANIZER

ENDO BLOCK



Graduated holes of different lengths provide an accurate, quick method of adjusting root canal instruments at the required length.

SPECIFICATIONS

Code	Colour	Packing
REB	Blue, Grey	250 pcs / box

SECURE DISPENSER / ORGANIZER

ENDO BOXES




These organizers are made of plastic and have removable, colour-coded silicon cups for Gutta Percha points & Paper point along with storage space for burs / files / reamers / spreaders. Autoclavable up to 135°C

SPECIFICATIONS

Code	Shape	Packing
REBS36	Square, 36 holes	1 pc / box, 120 boxes / case
REBC30	Circular, 30 holes	1 pc / box, 64 boxes / case
REBR60	Rectangle, 60 holes	1 pc / box, 96 boxes / case

RETRACTORS AND EXPANDERS

CHEEK RETRACTOR



This accessory is for retraction of lips and cheek during a dental treatment. Cheek Retractors are durable, flexible and convenient to use. The smooth surface makes the patient comfortable.

SPECIFICATIONS

Code	Type	Packing	Autoclave
RCRC S, M, L	C type: Large, Medium & Small	2 pcs / box, 200 boxes / case	up to 134°C
RCRM S, M, L	M type: Large, Medium, Small	2 pcs / box, 200 boxes / case	up to 134°C
RCRO L, S	O type: Large, & Small	2 pcs / box, 147 boxes / case	up to 134°C
RCRY L, S	Y type: Large & Small	10 pcs / bag, 100 boxes / case	Not applicable

RETRACTORS AND EXPANDERS

MOUTH PROP



The upper and lower teeth can be easily propped up in the mouth. These props are made from comfortable and elastic medical plastic material, with no glue and peculiar smell. These props are helpful to avoid patients' fatigue. The tongue guards block tongue movement to eliminate the chance of accidental injury. These props are resistance to high temperatures of up to 121°C, and can also withstand high-pressure disinfection.

SPECIFICATIONS

Product	Code	Size	Packing
Mouth Prop	RMPL	Large	50 pcs / bag, 20 bags / case
Mouth Prop	RMPM	Medium	50 pcs / bag, 20 bags / case
Mouth Prop	RMPS	Small	50 pcs / bag, 20 bags / case

PROTECTIVES & ABSORPTIVE

BOUFFANT CAP



ORO's bouffant caps provide the perfect fit to provide full hair coverage, and to stop cross-contamination by hair. These caps are highly breathable, light weight and antistatic, with an elegant look.

SPECIFICATIONS

Code	Type	Packing
RBCB18, RBCB21	Strip 18", 21" Blue	100 pcs / bag 40 bags / case for 21 50 bags / case for 18
RBCW18, RBCW21	Strip 18", 21" White	100 pcs / bag 40 bags / case for 21 50 bags / case for 18

PROTECTIVES & ABSORPTIVE

FACE MASK



ORO's 3-ply gentle, non-absorbent and non-irritating mask remains cool on the skin for ultimate ease and comfort. The adjustable nose piece helps to customize the fit, while the ear loop ties are soft, stretchable and latex free. Follow size and BFE as 17.5cm X 9.5cm, BFE ≤ 95%~99.05%.

SPECIFICATIONS

Code	Type	Packing
RFM3EB	3 Ply Ear Loop Blue	50 pcs / box, 40 boxes / case
RFM3EP	3 Ply Ear Loop Pink	50 pcs / box, 40 boxes / case
RFM3EW	3 Ply Ear Loop White	50 pcs / box, 40 boxes / case
RFM3TB	3 Ply Tie on Blue	50 pcs / box, 40 boxes / case
RFM4EBL	4 Ply Ear Loop Black	50 pcs / box, 40 boxes / case
RFM4EP	4 Ply Ear Loop Pink	50 pcs / box, 40 boxes / case
RFM4EW	4 Ply Ear Loop White	50 pcs / box, 40 boxes / case
RFM4EB	4 Ply Ear Loop Blue	50 pcs / box, 40 boxes / case

PROTECTIVES & ABSORPTIVE

FACE SHIELD



Anti-fog full face protection shield protects the face and eyes from flying debris and blood / water splashes.

- The foam shield is disposable and can be used over eyeglasses. The fit is very comfortable with elastic.
- The re-usable shield is washable, easy to clean, and light weight for better comfort. The shields take quick position and fits over eyeglasses.

SPECIFICATIONS

Code	Type	Packing
RFSD	Face Shield Disposable 33x32 cm	10 pcs / bag, 20 bags / case
RFSSF10S	With frame & 10 shields	1 set / box, 30 boxes / case

PROTECTIVES & ABSORPTIVE

IMPLANT KIT



Implant kits are designed to maximize efficiency and simplify surgical procedures. Single time use all in one kit, We offer a range of surgical kits specially configured to meet doctors' needs for virtually all clinical scenarios. The kit layout features easy-to-follow protocol flow, to assist in economizing time and effort for every phase of surgery.

SPECIFICATIONS

Code	Size	Packing
RIK	Instrument table cover 80cmx100cm	20 Kits/Case
	3 Ply Earloop face mask	
	Bouffant Cap	
	Tie back surgical gown	
	Shoes Cover	
	Towel 25cmx25cm	
	U non-woven drape with adhesive (140cmx170cm)	
	Non-woven sponges 3x3x4ply	
Tubing sleeve (plastic sheet with adhesive 8.5cmx120cm)		

PROTECTIVES & ABSORPTIVE

PATIENT BIB



ORO's 3-ply patient bib is made of soft tissue and poly of highly absorbent nature. The bonding of 2-ply tissues and 1-ply poly work securely together to eliminate separation of layers, and the reinforced edges provide added strength and durability. The tissues resist ripping even when wet, and can therefore double as a counter top drape for soaking up fluids.

Twisted clips are suitable to be used with the bibs while operating on the patient. The sufficient length ensures comfort.

SPECIFICATIONS

Product	Code	Colour	Packing
patient Bib	RPBB	Blue	125 pcs / bag, 4 bags / case
patient Bib	RPBG	Green	125 pcs / bag, 4 bags / case
patient Bib	RPBL	Lavender	125 pcs / bag, 4 bags / case
patient Bib	RPBO	Orange	125 pcs / bag, 4 bags / case
patient Bib	RPBP	Pink	125 pcs / bag, 4 bags / case
patient Bib	RPBW	White	20 pcs / bags, 40 bags / case

PROTECTIVES & ABSORPTIVE

DENTAL APRON



Dental apron help to protect skin and clothing from damage such as liquid splashing in the torso area. The apron has good absorbing capacity and its pattern is perfect to channelize liquid downward.

SPECIFICATIONS

Code	Type	Packing
RPANWTB	Non Woven Tie Back	100 pcs / bag, 10 bags / case
RPAVTB	Vinyl Tie Back	100 pcs / bag, 10 bags / case

PROTECTIVES & ABSORPTIVE

PROTECTIVE GLASS



Protective glasses are specially designed to protect the eyes from injury. These glasses are anti-fog, anti-hard light, and anti-ultraviolet ray. With a smooth surface they provide great comfort. The glasses also feature an anti-scratch material, which offers excellent protection and opal view.

SPECIFICATIONS

Code	Type	Packing
REWO	Orange	12 pcs / box, 12 boxes / case
REWC	Clear	12 pcs / box, 12 boxes / case

PROTECTIVES & ABSORPTIVE

SMS LAB COAT & JACKET



SMS protective lab coats and jackets come with pockets and knitted wrist cuffs. These are considered as elegant, professional wear, which also give better protection. The soft, breathable and fluid-resistant material offers matchless comfort.

SPECIFICATIONS

Code	Size	Colour	Packing
RJSBL RJSBM, RJSBS	S, M, L	Blue, Pink	10 pcs / bag, 100 bags / case
RLCSBL, RLCSBM, RLCSBS, RLCSBXL RLCSBXXL	XL, XXL L, M, S	Blue, Pink	10 pcs / bag, 100 bags / case

PROTECTIVES & ABSORPTIVE

SURGICAL GOWN



ORO's economical, disposable barrier gowns are made of fluid-repelling spun bond polypropylene. These gowns offer basic protection against light fluid spray, dirt, and dry particulates. This popular fabric is economical, breathable, light weight and flexible. The multi-purpose isolation gown features elasticized wrists, neck and waist ties, and overlapping back panels.

SPECIFICATIONS

Code	Description	Packing
RIGNS	Non Sterile	10 pcs / bag, 10 bags / case
RIGS	Sterile	1 pcs / bag, 8 bags / case

PROTECTIVES & ABSORPTIVE

COTTON ROLL AND DISPENSER



Cotton rolls are highly absorbent and made using 100% pure virgin white cotton with long fibres. These rolls are highly flexible, conform easily and maintain their shape during use. The absorption capacity is 10 times of its normal weight.

A single wrap cotton roll fits perfectly in the dispenser which gives it a smart look, while maintaining safety from cross-contamination. It has a rubber base for a sturdy non-slip surface grip with a transparent body.

SPECIFICATIONS

Product	Code	Size	Packing
Cotton Roll	RCR1X500	500 PCS	50 pcs / wrap, 10 wraps / bag, 80 bags / case
Cotton Roll	RCR2X500	500 PCS	50 pcs / wrap, 10 wraps / bag, 70 bags / case
Cotton Roll	RCR3X500	500 PCS	50 pcs / wrap, 10 wraps / bag, 60 bags / case
Cotton Roll	RCDP	Dispenser press type	1pcs / box, 40 boxes/ case

PROTECTIVES & ABSORPTIVE

SPONGES



ORO's non-sterilized non-woven sponges are soft, smooth and non-linting. The material offers maximum absorbency and little adhesion to wounds with pure white quality. Soft and comfortable to use either dry or moistened, these sponges absorb 10 times their weight.

SPECIFICATIONS

Code	Size	Packing
RSPGN4P2X2	4ply, 2"x2"	200 pcs / sleeve, 60 sleeves / case
RSPGN4P3X3	4ply, 3"x3"	200 pcs / sleeve, 24 sleeves / case
RSPGN4P4X4	4ply, 4"x4"	200 pcs / sleeve, 24 sleeves / case

SAFE KEEPERS

BUR HOLDER



Bur holders will organize your instruments for use during the procedure and while sterilization process. You can use these holders during procedures to give you quick and easy accessibility to your burs. Varying hole sizes can accommodate different types of shanks and cutting heads.

SPECIFICATIONS

Code	Shape & Shot	Packing
RBHAC60	Circular, 60 holes, Autoclavable	1 pc / box, 60 boxes / case
RBHAS24	Square, 24 holes, Autoclavable	1 pc / box, 100 boxes / case

SAFE KEEPERS

DENTURE BOX



ORO's denture box is a robust storage container for hygienic storage of dentures, braces, gum shields and other orthodontic appliances. The boxes are made of superior plastic to keep them healthy and clean.

SPECIFICATIONS

Code	Description	Colour	Packing
RDBB, RDBG, RDBLB, RDBR,	Denture Box	Blue & Green, Light Blue, Red	12 pcs / box, 16 boxes / case 296 pcs / case
RDBT	Denture Box with Tray	Dark Green, Green with White box	180 pcs / case

SAFE KEEPERS

RETAINER BOX



ORO's retainer box is made of sturdy plastic to protect the patient's retainer from bacteria or breakage. The retainer box includes vent holes or a hinged lid.

SPECIFICATIONS

Code	Colour	Packing
RRBS	Yellow, Clear Blue, Clear Green, Rose Red, White, Pink	12 pcs / box, 16 boxes / case or 450 pcs / case
RRBT	Yellow, Blue, Green, Rose Red, White	12 pcs / box, 16 boxes / case or 400 pcs / case

MISCELLANEOUS

GLASS DAPPEN DISH



These are dense glass vessels for mixing dental medicaments or fillings. The thick and heavy base makes the dish stable and offers protection from spill. It is ideal for mixing acrylics, amalgam, resins, sealants and other materials.

SPECIFICATIONS

Code	Colour	Packing
RDDG	Amber, Blue, Green, Clear	12pcs / box, 18 boxes / case

MISCELLANEOUS

PRE-BENT TIPS



Pre-bent tips are disposable, blunt-end needle tips designed to dispense etchants, resins and flowable composites cleanly and easily. These tips are offered in several different gauge sizes to accommodate most materials and procedures.

SPECIFICATIONS

Code	Gauge	Packing
RPBT18	18ga (cavity liners, bases and composites)	100 pcs / bag
RPBT19	19ga (flowable composites and sealants)	100 pcs / bag
RPBT20	20ga (flowable composites)	100 pcs / bag
RPBT21	22ga (sealants and cements)	100 pcs / bag
RPBT22	25ga (etchants)	100 pcs / bag

MISCELLANEOUS

WOODEN WEDGES



Eco-friendly wooden wedges are anatomically shaped for better and safer procedures. We offer 6 wedges which can be trimmed without splintering, for a perfect custom fit. The wedges' colour-coding agent does not recolour" the tooth or restoration upon removal.

SPECIFICATIONS

Code	Size (mm)	Colour	Packing
RWWA	11mm x 3.5mm x 1.5mm x 1.3mm	Assorted	100 pcs / box, 36 boxes / case
RWWB	14mm x 4mm x 1.5mm x 1.5mm	Blue	100 pcs / box, 36 boxes / case
RWWG	14mm x 2mm x 1.5mm x 1.5mm	Green	100 pcs / box, 36 boxes / case
RWWO	16mm x 4mm x 1.5mm x 1.5mm	Orange	100 pcs / box, 36 boxes / case
RWWP	14mm x 4mm x 2mm x 1.5mm	Pink	100 pcs / box, 36 boxes / case
RWWW	16mm x 4.5mm x 2mm x 1.5mm	White	100 pcs / box, 36 boxes / case
RWWY	16mm x 4.5mm x	Yellow	100 pcs / box, 36 boxes / case

MISCELLANEOUS

PROPHY BRUSH & CUP



ORO's prophy brush reaches into recesses where a normal brush can't reach. The variant shape of the brush removes heavy stains from occlusal surfaces and hard-to-reach inter-proximal spaces. Prepares surface for sealant and orthodontic procedure. It is made up of nylon, and is latex-free.

SPECIFICATIONS

Code	Type	Image	Packing
RPB-330	latch type flat white		100 pcs / box, 100 boxes / case
RBP-340	latch type tapered white		100 pcs / box, 100 boxes / case
RBP-350	latch type flat black		100 pcs / box, 100 boxes / case
RBP-360	latch type tapered black		100 pcs / box, 100 boxes / case
RBP-370	latch type round white		100 pcs / box, 100 boxes / case
RPB-330A	latch type flat assorted		100 pcs / box, 100 boxes / case
RPB-370A	latch type round assorted		100 pcs / box, 100 boxes / case
RPC	cup latch type style 6 rib		

MISCELLANEOUS

X-RAY FILM HOLDER



X-ray holder takes quick and precise anterior and posterior bisecting angle radiographs, and eliminates manual retention on the film packet by patient. It enhances patient comfort with smooth edges.

SPECIFICATIONS

Code	Colour	Packing
RXFH	Beige	4 pcs / box, 200 boxes / case



REACH GLOBAL INDIA PVT. LTD.

Head Office:

Office No. 1 & 2, Survey No. 64/1, Wing-A, First Floor, Sai Raj Regency,
Pimple Saudagar, Pune, State - Maharashtra (India), Pin Code - 411027

Manufacturing Unit:

Plot No I-22, Gat No. 1064, Chordia Industrial Park, Naigaon, Shirwal,
Taluka-Khandala, Dist. Satara, State-Maharashtra (India) Pin Code - 412801

info@orohealthcare.com
www.orohealthcare.com