



# HEALTHCARE RANGE



**MADE IN INDIA**



[www.orohealthcare.com](http://www.orohealthcare.com)



**Committed to build safe and better world for healthcare industry**  
Think sterilization and safety | Think ORO



# REACH GLOBAL INDIA PVT. LTD.

Redefining standards in the Healthcare Industry

Reach Global India Pvt. Ltd. is a reputed name in the Healthcare industry for product quality, availability and technical support with best-in-class service. Over the years, Reach Global has consistently kept itself abreast of the latest technologies and stayed ahead of the curve in the industry.

We believe in offering one-stop solution to the customer. Reach Global India Pvt. Ltd. always do their best to achieve 'value for the money' products, services by maintaining operations cost effective and at the same time by keeping in mind customer requirements towards increasing compliance by keeping 'customer at centre'. Through our well spread of extended arm network across PAN India and Global, we try to remain closer to customer.

With a commitment to deliver 'Your Infection Guard', Reach Global India Pvt. Ltd. is pillared on the values of high quality, transparent dealings and flawless distribution to stay at the forefront of meeting each and every requirement of the Healthcare industry in India and overseas.

## CAPABILITIES

- CE Certified Sterilization packaging products
- EN ISO 13485 certified and IMDR compliant state-of-the-art manufacturing facility
- ISO- 8 Certified cleanroom environment
- USFDA approved safety Products
- In-house production & In-stock product availability
- Dedicated warehouse management system





## VISION

While retaining our current reputation of "your infection guard" Our vision is to become 'leading manufacturer and supplier' of barrier packaging & and associated products for controlled area.

## MISSION

Our Sterilization packaging products, Protective clothing and supplementary products are targeted to cater the needs of 'Healthcare segment' for their CSSD, Production, Quality Control & Assurance sections. We aim to support their needs by offering suitable products manufactured at 'state-of-the art' facility with desire to maintain the 'best availability'




## OUR VALUES

### BELIEF

Our strong belief is our people and we work by taking ownership to deliver our promises persistently

### FUTURISTIC

We plan ourselves by targeting long term and sustainable success



### UNIFY

We encourage and challenge each other with respect to achieve our common objective with fun and by helping each other unified

### UNDERSTANDING

We are open and have courage to try new and learn from our mistakes. we get happiness in developing partnership with our customer by rising understanding of small things to offer a solution



## STERILIZATION PACKAGING **Page**

Sterilization Reel Flat .....	6
Sterilization Reel Gussetted.....	7
Sterilization Reel HDPE.....	8
Sterilization Reel Plasma.....	9
Sterilization Pouch Self Seal & Heat Seal .....	10
Wrap Sheets SMS .....	11
& Wrap Sheets Crepe	

## STERILITY MONITORING **Page**

Bowie Dick Test Pack .....	12
& Helix Strip Indicator	
Chemical Indicator Type 4/5/6.....	13
Indicator Tapes & Process Indicator Label.....	14

## SAFETY WEAR **Page**

Shoe Cover & SMS Jacket .....	15
SMS Labcoat and Isolation Gown, Sterile Surgeon Gown & Reinforce Gown.....	16
Surgeon Cap & Bouffant Cap.....	17
Face Mask.....	18
Disposable Coverall & Reusable Coverall.....	19

# STERILIZATION PACKAGING



## Medical Grade Reel Flat

Medical grade sterilization reels allow sterilization of packaging contents while maintaining sterility until the package is opened. Ensuring storage, handling and transportation delivery to avoid any cross contamination

- Manufactured using virgin medical grade imported paper in compliant with ISO111607-1
- Triple line sealing makes it easy to peel
- Wide range of standard sizes with flexibility to offer custom made sizes
- Tinted transparent, reinforced multiayer co-polymer PET/CPP film
- Excellent steam permeability and microbial barrier properties
- Recomend for light weight products not having sharp edges

*\*Sharp object must have covers to avoid punctures*



CODE	PRODUCT NAME	PKG (Reels/Ctn)
RRM50	STERILIZATION REEL 50x200M	10
RRM55	STERILIZATION REEL 55x200M	10
RRM75	STERILIZATION REEL 75x200M	8
RRM100	STERILIZATION REEL 100x200M	6
RRM125	STERILIZATION REEL 125x200M	8
RRM150	STERILIZATION REEL 150x200M	4
RRM200	STERILIZATION REEL 200x200M	4
RRM250	STERILIZATION REEL 250x200M	4
RRM300	STERILIZATION REEL 300x200M	2
RRM350	STERILIZATION REEL 350x200M	2
RRM400	STERILIZATION REEL 400x200M	2
RRM450	STERILIZATION REEL 450x200M	2
RRM500	STERILIZATION REEL 500x200M	2



*\*Customized sizes available on demand*



# STERILIZATION PACKAGING




## Gusseted Reel

Medical grade sterilization reels allow sterilization of packaging contents while maintaining sterility until the package is opened. Ensuring storage, handling and transportation delivery to avoid any cross contamination

- Manufactured using virgin medical grade imported Paper in compliant with ISO111607-1
- Triple line sealing makes it easy to peel
- Wide range of standard sizes with flexibility to offer custom made sizes
- Tinted transparent, reinforced multiayer co-polymer PET/CPP film
- Excellent steam permeability and microbial barrier properties
- Recomendaded for light weight products not having sharp edges

*\*Sharp object must have covers to avoid punctures*



CODE	PRODUCT NAME	PKG (Reels/Ctn)
RRMG50	STERILIZATION REEL GUSSETED 50x100M	10
RRMG55	STERILIZATION REEL GUSSETED 55x100M	10
RRMG75	STERILIZATION REEL GUSSETED 75x100M	8
RRMG100	STERILIZATION REEL GUSSETED 100x100M	6
RRMG125	STERILIZATION REEL GUSSETED 125x100M	8
RRMG150	STERILIZATION REEL GUSSETED 150x100M	4
RRMG200	STERILIZATION REEL GUSSETED 200x100M	4
RRMG250	STERILIZATION REEL GUSSETED 250x100M	4
RRMG300	STERILIZATION REEL GUSSETED 300x100M	2
RRMG350	STERILIZATION REEL GUSSETED 350x100M	2
RRMG400	STERILIZATION REEL GUSSETED 400x100M	2
RRMG450	STERILIZATION REEL GUSSETED 450x100M	2
RRMG500	STERILIZATION REEL GUSSETED 500x100M	2

EN868-5

**\*Customized sizes available on demand**

# STERILIZATION PACKAGING



## HDPE Reel

Nonwoven-Spunbond HDPE widely preferred for its outstanding strength, durability puncture & tear resistance

- Excellent bacterial barrier properties against bacteria, micro organism their spores
- Puncture resistant and high tear resistant
- Excellent seal and aseptic presentation
- Clear and accurate indicator color transition
- Safe to use for heavy weight articles
- Provides less particle generation while opening the package, preventing contamination with the fibers
- Recommended for Steam & ETO



CODE	PRODUCT NAME	PKG (Reels/Ctn)
RRHD50	ORO STERILIZATION REEL HDPE 50x200M	10
RRHD55	ORO STERILIZATION REEL HDPE 55x200M	10
RRHD75	ORO STERILIZATION REEL HDPE 75x200M	8
RRHD100	ORO STERILIZATION REEL HDPE 100x200M	2
RRHD125	ORO STERILIZATION REEL HDPE 125x200M	2
RRHD150	ORO STERILIZATION REEL HDPE 150x200M	2
RRHD200	ORO STERILIZATION REEL HDPE 200x200M	1
RRHD250	ORO STERILIZATION REEL HDPE 250x200M	1
RRHD300	ORO STERILIZATION REEL HDPE 300x200M	1
RRHD350	ORO STERILIZATION REEL HDPE 350x200M	1
RRHD400	ORO STERILIZATION REEL HDPE 400x200M	1
RRHD450	ORO STERILIZATION REEL HDPE 450x200M	1
RRHD500	ORO STERILIZATION REEL HDPE 500x200M	1

EN868-5

\*Customized sizes available on demand



# STERILIZATION PACKAGING



## Plasma Reel

Plasma reels are made of nonwoven HDPE fabric which provides outstanding strength, durability, puncture & tear resistance

- Excellent bacterial barrier properties against bacteria, micro organism & their spores
- Extremely high resistance to penetration and fibre-free peeling
- Excellent seal and aseptic presentation
- Clear and accurate indicator color transition
- Long shelf life which can maintain sterility for at least five years, if package integrity is not compromised



CODE	PRODUCT NAME	PKG (Reels/Ctn)
RRP5070	STERILIZATION REEL PLASMA 50x70M	10
RRP5570	STERILIZATION REEL PLASMA 55x70M	10
RRP7570	STERILIZATION REEL PLASMA 75x70M	8
RRP10070	STERILIZATION REEL PLASMA 100x70M	6
RRP12570	STERILIZATION REEL PLASMA 125x70M	8
RRP15070	STERILIZATION REEL PLASMA 150x70M	4
RRP20070	STERILIZATION REEL PLASMA 200x70M	4
RRP25070	STERILIZATION REEL PLASMA 250x70M	4
RRP30070	STERILIZATION REEL PLASMA 300x70M	2
RRP35070	STERILIZATION REEL PLASMA 350x70M	2
RRP40070	STERILIZATION REEL PLASMA 400x70M	2
RRP45070	STERILIZATION REEL PLASMA 450x70M	2
RRP50070	STERILIZATION REEL PLASMA 500x70M	2

EN868-5

\*Customized sizes available on demand

# STERILIZATION PACKAGING



## Sterilization Pouches Heat-Seal & Self-Seal

ORO sterilization pouches are specially designed for convenient and safe packaging of items to be sterilized, ensuring sterility of the articles until package integrity is not compromise. The corner tack seals prevents corners from curling after sterilization. Available in self-seal and heat-seal.

- Chevron triple seal provides easy peel experience, while withstanding the sealing strength
- Available in Medical grade for steam & EO application
- Available in Hyapak & HDPE for plasma & EO application
- Packaging: 200pcs/Bag



EN868-5

## Sterilization Self-Seal Pouch Sizes

70 X 255 MM	75 X 280 MM	90 X 165 MM	100 X 280 MM	190 X 363 MM
200 X 280 MM	330 X 450 MM	200 X 345 MM	250 X 345 MM	250 X 400 MM
300 X 400 MM	400 X 500 MM	500 X 500 MM	90 X 260 MM	134 X 283 MM

## Sterilization Heat-Seal Pouch Sizes

70 X 230 MM	75 X 254 MM	90 X 138 MM	100 X 253 MM	190 X 337 MM
200 X 254 MM	75 X 240 MM	200 X 315 MM	250 X 376 MM	300 X 372 MM
330 X 426 MM	400 X 475 MM	500 X 473 MM	90 X 233 MM	134 X 253 MM

**\*Customized sizes available on demand**

*\* For product code and packaging please refer price list or contact ORO sales representative*

# STERILIZATION PACKAGING




## Wrap Sheets SMS

SMS are made with layers of spunbound & meltblown overlapping each other

- Single use sheets, recommended for heavy and bulky products
- Excellent water-repellent material
- Smart Packaging for convenient usage
- Offers high tensile strength
- In compliant with ISO11607-1

**Customized sizes available in 35 GSM, 43 GSM, 60 GSM etc**




## Wrap Sheets Crepe

Medical grade crepe paper enables high resistance and excellent drapability for articles to be sterilized. It is water repellent and works as the perfect barrier against contamination.

- Ecofriendly widely used & accepted packaging material.
- Soft fabrics provide most secure packaging
- High-quality light weight material
- Smart packaging for convenient usage
- Wrap sheets crepe are processed from the direct virgin pulp of tree
- In compliant with ISO11607-1
- Recommended for light weight products not having sharp edges

**Customized sizes available in 45 GSM, 60 GSM etc**



\*Colours available

40CM*40CM	50CM*50CM	60CM*60CM	75CM*75CM
80CM*80CM	90CM*90CM	100CM*100CM	120CM*120CM

\* For product code and packaging please refer price list or contact ORO sales representative




# STERILITY MONITORING



## Bowie Dick Test Pack

Bowie & Dick test Pack for air-removal test are manufactured as per the ISO 11140-1 and compliant with ISO11140-4.

- Used to detect air leakages & residual condensate in Autoclave machine
- Precise and sharp indicator color transition from blue to pink change provides distinctive pass/fail result
- Each test pack has Type-1 indicator dot on top for easy identification of processed or unprocessed packages
- 100 % lead free and heavy metal free indicator ink



PRODUCT CODE	PRODUCT NAME	PACKAGING
RSMBD-121/8.3	ORO BOWIE DICK TEST PACK TYPE-2 121°C-8.3MIN/134°C-3.5MIN	120 Pcs/ctn
RSMBD-121/10	ORO BOWIE DICK TEST PACK TYPE-2 121°C-10MIN/134°C-3.5MIN	
RSMBD-121/11	ORO BOWIE DICK TEST PACK TYPE-2 121°C-11MIN/134°C-3.5MIN	
RSMBD-121/15	ORO BOWIE DICK TEST PACK TYPE-2 121°C-15MIN/134°C-3.5MIN	
RSMBDS	ORO BOWIE DICK TEST SHEET-A4	50 Pcs/Bag

## Helix Test Strips

- Ensures air removal and steam penetration in pre-vacuum sterilizers.
- EN ISO 11140-1:2014 and EN ISO 11140-4
- Precise and sharp process indicator colour change pink to dark brown/black indicates successful sterilization

PRODUCT NAME	PRODUCT CODE	PACKAGING
Helix Strip -4	RBDHS4S	10 Strips/sheet
Helix Strip -6	RBDHS6S	25 Sheets/Box



# STERILITY MONITORING




## Chemical Indicator Strips

Chemical Indicator strips are essential tools for monitoring sterilization. They provide quick, reliable visual confirmation for the critical parameters like temperature, pressure, or chemical exposure have been met. ISO 11140-1 compliant

- Self adhesive strips
- For Steam Application available in Type 4,5&6
- For EO application available in Type 5 which ensures accurate monitoring in EO sterilizers
- For Plasma application available in Type 4, confirms effective plasma sterilization conditions
- 100% lead free and heavy metal free indicator ink

PRODUCT NAME	PRODUCT CODE
ORO CHEMICAL INDICATOR STEAM TYPE-4 134°C-3.5MIN/121°C-15min	RSMCIST4-121/15
ORO CHEMICAL INDICATOR STEAM TYPE-4 134°C-7MIN/121°C-20min	RSMCIST4-121/20
ORO CHEMICAL INDICATOR STEAM TYPE-5 134°C-3.5MIN/121°C-15MIN	RSMCIST5-121/15
ORO CHEMICAL INDICATOR STEAM TYPE-5 134°C-7MIN/121°C-20MIN	RSMCIST5-121/20
ORO CHEMICAL INDICATOR STEAM TYPE-6 134°C- 3.5MIN/121°C-15MIN	RSMCIST6-121/15
ORO CHEMICAL INDICATOR STEAM TYPE-6 134°C-7MIN/121°C-20MIN	RSMCIST6-121/20
ORO CHEMICAL INDICATOR EO TYPE-5	RSMCIET5
ORO CHEMICAL INDICATOR PLASMA TYPE-4	RSMCIPT4
PACKAGING	200 Strips/Box





# STERILITY MONITORING



## Indicator Tapes

ORO provides the accurate solution for safe sterility monitoring of sterilization wraps using these indicator tapes. These are non-hazardous tapes available for all kind steam, EO and plasma sterilization methods. The sharp indicator inks give a precise color change when meeting the prescribed sterilization method. The strong adhesive provides a firm & secure sealing keeping the integrity of the package intact, giving a comfortable & clean peeling leaving no residues behind while unpacking the package in aseptic area.

- High quality non-toxic adhesive
- Compatible to use over SMS wrap sheets, crepe sheets paper & film
- Available in 19mm & 25mm width with 50m length



Code	RSMSIT	RSMEIT	RSMPIT
Type	Steam	EO	Plasma
PKG	48 Pcs/Ctn		

## Process Indicator Label

ORO 3-line type-1 process indicator label on exposure to the sterilization process, helping to differentiate between processed and unprocessed loads

- 100% lead free and heavy metal free indicator ink
- Precise & sharp process indicator colour change
- High quality non toxic adhesive
- Packaging 1 roll x 520 nos. labels
- ISO 11140-1 compliant



Code	RPILS	RPILE	RPILP
Type	Steam	EO	PLASMA
PKG	60 Rolls/Ctn		

*\*For Label Gun please contact ORO sales representative*



# SAFETY WEAR



## Shoe Cover

- The ORO universal size Shoe Cover made of Non-Woven or LDPE Fabric.
- It is water resistant, anti-static & Anti-Skid.
- The Shoe cover gives great comfort for long hours working with great comfort & extreme hygiene.



CODE	PRODUCT NAME	PKG	
RSHCN	SHOE COVER NON-WOVEN	100 Pcs / Bag	28 Bag / Ctn
RSHCL	SHOE COVER LDPE	100 Pcs / Bag	50 Bag / Ctn

## SMS Jacket

SMS Jackets are lightweight, breathable protective garments made from Spunbond-Meltblown-Spunbond material. They provide excellent barrier protection against fluids and contaminants while ensuring comfort and durability, making them ideal for healthcare use.

- Knitted cuff & caw for comfort and complete enclosure
- Soft breathable & fluid resistant fabric assuring all day protection
- Pockets & snap closure



CODE	SIZES	PKG
RJSB	SMALL, MEDIUM, LARGE, XL	10 Pcs / Bag, 10 Bags/Ctn
RJSL	SMALL, MEDIUM, LARGE, XL	
RJSP	SMALL, MEDIUM, LARGE, XL	

# SAFETY WEAR



## Surgeon Cap

The disposable surgical caps are meant to be used by doctors and nurses while conducting surgeries. ORO caps are manufactured using optimum quality fabrics, which gives great comfort & highly breathable. Soft fabric with smooth texture makes the caps light weighted Backside elastic for easy size adjustment.



CODE	PRODUCT NAME	TYPE	PACKAGING
RSCNW	SURGEON CAP	BLUE	100 Pcs / Bag

## Bouffant Cap

ORO's bouffant caps provide the perfect fit to provide full hair coverage, and to stop cross-contamination by hair. These caps are highly breathable, light weight and antistatic, with an elegant look.



CODE	PRODUCT NAME	TYPE	PACKAGING
RBCS18B, RBCS21B	BOUFFANT CAP BLUE	STRIP 18" & 21"	100 Pcs / Bag
RBCS18W, RBCS21W	BOUFFANT CAP WHITE	STRIP 18" & 21" WHITE	100 Pcs / Bag

# SAFETY WEAR



## SMS Labcoat

ORO's Lab Coats are considered as elegant, professional wear, which also give better protection. The soft, breathable and fluid-resistant material offers matchless comfort.

- Made of non-woven SMS
- Knitted cuff & collar for comfort and complete enclosure
- Soft breathable & fluid resistant fabric assuring all day protection
- Pockets & snap closure
- Available sizes- S, M, L, XL, XXL

CODE	SIZE	PACKAGING
RLCSB	S, M, L, XL, XXL	100 Pcs / Bag

**\*Customized sizes available on demand**

## Isolation Gown, Sterile Surgeon Gown & Reinforce Gown

The non-woven gowns are highly breathable providing best protection

- Excellent fluid resistant nonwoven fabric
- High microbial barrier
- Great breathability giving all day comfort
- Ultimate protection against dirt & dry particulates
- Ideal for sterile area usage

CODE	RIGR	RIGS	PKG	RIG	PKG
TYPE	REINFORCED STERILE	STERILE SURGEON GOWN	1 Pcs/Bag	ISOLATION GOWN	1 Pcs/Bag



*\*Please check for MOQ  
of colors except blue*





# SAFETY WEAR



## Face Mask

ORO non-absorbent & non-irritating masks remains cool on the skin providing ultimate comfort

- Fluid Resistant mask
- Adjustable nose strip
- Soft ear loop for comfortable wear
- TYPE II IR mask confirming to EN14683:2019+AC:2019
- $\geq 98.4\%$  BFE, Superior Melt blown
- Ear loop & Tie-on variants available
- Available in 3PLY & 4PLY



\*Colour available for 3 Ply & 4 Ply



CODE	PRODUCT NAME	PACKAGING
RFM3E	FACE MASK 3 PLY EAR LOOP (BLUE, BLACK, PINK, WHITE)	50 Pcs / box
RFM4E	FACE MASK 4 PLY EAR LOOP (BLUE, BLACK PINK, WHITE)	50 Pcs / box
RFM3T	FACE MASK 3 PLY TIE ON (BLUE, BLACK PINK, WHITE)	50 Pcs / box
RFM4T	FACE MASK 3 PLY TIE ON (BLUE, BLACK PINK, WHITE)	50 Pcs / box

# SAFETY WEAR



## Disposable Coverall

Disposable Coverall is designed to cover the entire body and avoid cross-contaminations in the controlled environment and process carried out in the controlled. Single use disposable coverall act as a barrier for controlled contamination by human activities.

- Suitable for controlled environment
- Pre-sterile for quick and easy access
- Made of breathable PP+PE Non-Woven Fabric
- EN13795-2 complaint

*\*For details contact ORO sales representative*



## Reusable Coverall

Reusable Coverall plays key role in the contamination control of clean room environment. Human is the sole source of contamination in clean room environments apart from the material. Human skin carries various types of microorganisms like Bacteria, Fungal spores and viruses which can be carried with the scale flax and contaminate the clean room environment. ORO Reusable Coverall are manufactured by using 98% Polyester and 2% carbon Fibers, which have antistatic properties. Reusable Coverall are designed and packed in such a way that aseptic operations can be fulfilled. 3 Piece Reusable coverall covers complete body from top to toe and avoid the contaminations through particles shading in the cleanroom. ORO Reusable Coverall are breathable to achieve comfortable operations.

*\*For details contact ORO sales representative*





## Reach Global India Pvt. Ltd

- Reach Global India Private Limited.  
Office No. 1&2, Survey No. 64/1, Wing-A, First Floor,  
Sai Raj Regency, Pimple Saudagar, Pune: 411027
- Plot No 1-22, Gat No. 1064, Chordia Industrial Park,  
Naigaon, Shirwal, Taluka-Khandala,  
Dist. Satara, State-Maharashtra (India)  
Pin Code - 412801

Customer Care: +91 9370 788 850/ +91-(020) 27400987  
info@orohealthcare.com

