



HEALTHCARE RANGE

MEDICAL DEVICES



MADE IN INDIA

 www.orohealthcare.com





Committed to build safe and better world for healthcare industry
Think sterilization and safety | Think ORO



REACH GLOBAL INDIA PVT. LTD.

Redefining standards in the Healthcare Industry

Reach Global India Pvt. Ltd. is a highly reputed name in the Healthcare industry in terms of product quality and financial performance with best-in-class service. Over the years, Reach Global has consistently kept itself abreast of the latest technologies and stayed ahead of the curve in the industry.

Operating as a true one-stop solution, Reach Global has always provided its customer network a basket of high-range performance products and solutions at an optimum cost. Besides demonstrated performance, durability, low maintenance and safe operations, the brand has earned itself the niche of also being rated high on cost-efficiency. Thereby delivering optimum value with thorough professionalism in business operations.

With a commitment to deliver 'Your Infection Guard', Reach Global India Pvt. Ltd. - is pillared on the values of high quality, transparent dealings and flawless distribution to stay at the forefront of meeting each and every requirement of the dental industry in India and overseas.

CAPABILITIES

- CE Certified Sterilization packaging product
- EN ISO 13485 certified and IMDR compliant state-of-the-art manufacturing facility
- ISO- 8 Certified cleanroom environment
- USFDA approved safety Products
- In-house production & In-stock product availability
- Dedicated warehouse management system



VISION

While retaining our current reputation of "your infection guard" Our vision is to become 'leading manufacturer and supplier' of barrier packaging & and associated products for controlled area.



MISSION

Our Sterilization packaging products, Protective clothing and supplementary products are targeted to cater the needs of 'Healthcare segment' for their CSSD, Production, Quality Control & Assurance sections. We aim to support their needs by offering suitable products manufactured at 'state-of-the art' facility with desire to maintain the 'best availability'

OUR VALUES



BELIEF

Our strong belief is our people and we work by taking ownership to deliver our promises persistently



UNIFY

We encourage and challenge each other with respect to achieve our common objective with fun and by helping each other unified



FUTURISTIC

We plan ourselves by targeting long term and sustainable success



UNDERSTANDING

We are open and have courage to try new and learn from our mistakes. we get happiness in developing partnership with our customer by rising understanding of small things to offer a solution.



STERILIZATION PACKAGING

Page

Sterilization Reel Flat.....	6
Sterilization Reel Gussetted.....	7
Sterilization Reel HDPE.....	8
Sterilization Hypak Reel	9
Sterilization Pouch Self Seal & Heat Seal.....	10
Cut Sheets & Rolls Sheet	11

STERILITY MONITORING

Page

Bowie Dick Test Pack.....	12
& Helix Strip Indicator	
Chemical Indicator Type 4/5/6.....	13
Indicator Tapes & Process Indicator Label.....	14

SAFETY WEAR

Page

Shoe Cover & SMS Jacket	15
SMS Labcoat and Isolation Gown, Sterile Surgeon Gown & Reinforce Gown.....	16
Surgeon Cap & Bouffant Cap.....	17
Face Mask.....	18
Disposable Coverall & Reusable Coverall.....	19

STERILIZATION PACKAGING



MEDICAL GRADE REEL FLAT

Medical grade sterilization reels allow sterilization of packaging contents while maintaining sterility until the package is opened. Ensuring storage, handling and transportation delivery to avoid any cross contamination

- Manufactured using virgin medical grade imported paper in compliant with ISO111607-1
- Triple line sealing makes it easy to peel
- Wide range of standard sizes with flexibility to offer custom made sizes
- Tinted transparent, reinforced multiayer co-polymer PET/CPP film
- Excellent steam permeability and microbial barrier properties
- Recommended for light weight products not having sharp edges

**Sharp object must have covers to avoid punctures*



CODE	PRODUCT NAME	PKG (Reels / Ctn)
RRM50	STERILIZATION REEL 50MMX200M	10
RRM55	STERILIZATION REEL 55MMX200M	10
RRM75	STERILIZATION REEL 75MMX200M	8
RRM100	STERILIZATION REEL 100MMX200M	6
RRM125	STERILIZATION REEL 125MMX200M	8
RRM150	STERILIZATION REEL 150MMX200M	4
RRM200	STERILIZATION REEL 200MMX200M	4
RRM250	STERILIZATION REEL 250MMX200M	4
RRM300	STERILIZATION REEL 300MMX200M	2
RRM350	STERILIZATION REEL 350MMX200M	2
RRM400	STERILIZATION REEL 400MMX200M	2
RRM450	STERILIZATION REEL 450MMX200M	2
RRM500	STERILIZATION REEL 500MMX200M	2

EN8685

*Customized sizes available on demand

STERILIZATION PACKAGING



GUSSETED REEL

Medical grade sterilization reels allow sterilization of packaging contents while maintaining sterility until the package is opened. Ensuring storage, handling and transportation delivery to avoid any cross contamination

- Manufactured using virgin medical grade imported paper in compliant with ISO111607-1
- Triple line sealing makes it easy to peel
- Wide range of standard sizes with flexibility to offer custom made sizes
- Tinted transparent, reinforced multilayer co-polymer PET/CPP film
- Excellent steam permeability and microbial barrier properties
- Recommended for light weight products not having sharp edges

**Sharp object must have covers to avoid punctures*



CODE	PRODUCT NAME	PKG (Reels / Ctn)
RRMG50	STERILIZATION REEL GUSSETED 50MMX200M	10
RRMG55	STERILIZATION REEL GUSSETED 55MMX200M	10
RRMG75	STERILIZATION REEL GUSSETED 75MMX200M	8
RRMG100	STERILIZATION REEL GUSSETED 100MMX200M	6
RRMG125	STERILIZATION REEL GUSSETED 125MMX200M	8
RRMG150	STERILIZATION REEL GUSSETED 150MMX200M	4
RRMG200	STERILIZATION REEL GUSSETED 200MMX200M	4
RRMG250	STERILIZATION REEL GUSSETED 250MMX200M	4
RRMG300	STERILIZATION REEL GUSSETED 300MMX200M	2
RRMG350	STERILIZATION REEL GUSSETED 350MMX200M	2
RRMG400	STERILIZATION REEL GUSSETED 400MMX200M	2
RRMG450	STERILIZATION REEL GUSSETED 450MMX200M	2
RRMG500	STERILIZATION REEL GUSSETED 500MMX200M	2

EN8685

*Customized sizes available on demand

STERILIZATION PACKAGING



HDPE REEL

Nonwoven-Spunbond HDPE widely preferred for its outstanding strength, durability, puncture & tear resistance

- Excellent bacterial barrier properties against bacteria, micro organism & their spores
- Puncture resistant and high tear resistant
- Excellent seal and aseptic presentation
- Clear and accurate indicator color transition
- Safe to use for heavy weight articles
- Provides less particle generation while opening the package, preventing contamination with the fibers
- Recommended for Steam & ETO



CODE	PRODUCT NAME	PKG
RRHD50	STERILIZATION REEL HDPE 50MMx200M	10
RRHD55	STERILIZATION REEL HDPE 55MMx200M	10
RRHD75	STERILIZATION REEL HDPE 75MMx200M	8
RRHD100	STERILIZATION REEL HDPE 100MMx200M	2
RRHD125	STERILIZATION REEL HDPE 125MMx200M	2
RRHD150	STERILIZATION REEL HDPE 150MMx200M	2
RRHD200	STERILIZATION REEL HDPE 200MMx200M	1
RRHD250	STERILIZATION REEL HDPE 250MMx200M	1
RRHD300	STERILIZATION REEL HDPE 300MMx200M	1
RRHD350	STERILIZATION REEL HDPE 350MMx200M	1
RRHD400	STERILIZATION REEL HDPE 400MMx200M	1
RRHD450	STERILIZATION REEL HDPE 450MMx200M	1
RRHD500	STERILIZATION REEL HDPE 500MMx200M	1

EN8685

*Customized sizes available on demand

STERILIZATION PACKAGING



Hypak® Reel

Hypak® — a metamaterial made from unique flash-spun technology. It's made from Nonwoven-Spunbond HDPE material.

- Allows exchange of air & gaseous sterilant during sterilization & acts as a barrier after sterilization
- Act as a selective barrier membrane, do not allow any virus, bacteria & microbe to pass through
- Polymeric composition is non-toxic and harmless, does not degrade upon exposure to various sterilization gases or high temperatures
- Provides minimum particle generation while opening healthcare devices, preventing contamination with the fibers
- Hypak material is in compliant with ISO11607-1



CODE	PRODUCT NAME	PKG
RRH200	ORO STERILIZATION REEL HYPAK 1075B 200 MM X 200M	4
RRH250	ORO STERILIZATION REEL HYPAK 1075B 250 MM X 200M	4
RRH300	ORO STERILIZATION REEL HYPAK 1075B 300 MM X 200M	2
RRH350	ORO STERILIZATION REEL HYPAK 1075B 350 MM X 200M	2
RRH400	ORO STERILIZATION REEL HYPAK 1075B 400 MM X 200M	2
RRH450	ORO STERILIZATION REEL HYPAK 1075B 450 MM X 200M	2
RRH500	ORO STERILIZATION REEL HYPAK 1075B 500 MM X 200M	2
RRH550	ORO STERILIZATION REEL HYPAK 1075B 550 MM X 200M	2
RRH600	ORO STERILIZATION REEL HYPAK 1075B 600 MM X 200M	2

EN8685

*Customized sizes available on demand

STERILIZATION PACKAGING



STERILIZATION POUCHES HEAT-SEAL & SELF-SEAL

ORO sterilization pouches are specially designed for convenient and safe packaging of items to be sterilized, ensuring sterility of the articles until package integrity is not compromise. The corner tack seals prevents corners from curling after sterilization. Available in self-seal and heat-seal.

- Chevron triple seal & thumb-notch provides easy peel experience, while with standing the sealing strength as well.
- Available in Medical Grade, Hypak and HDPE
- Packaging: 200pcs/Bag

EN8685



Sterilization Self-Seal Pouch Sizes

70 x 255 MM	75 X 280 MM	90 X 165 MM	100 X 280 MM	190 X 363 MM
200 X 280 MM	330 X 450 MM	200 X 345 MM	250 X 345 MM	250 X 400 MM
300 X 400 MM	400 X 500 MM	500 X 500 MM	90 X 260 MM	134 X 283 MM

Sterilization Heat-Seal Pouch Sizes

70 x 230 MM	75 X 254 MM	90 X 138 MM	100 X 253 MM	190 X 337 MM
200 X 254 MM	75 X 240 MM	200 X 315 MM	250 X 376 MM	300 X 372 MM
330 X 426 MM	400 X 475 MM	500 X 473 MM	90 X 233 MM	134 X 253 MM

*Customized sizes available on demand

* For product code and packaging please refer price list or contact ORO sales representative

HDPE SHEETS & ROLLS FOR BLISTER PACKAGING

We offer HDPE Sheets & Rolls specifically designed to meet the stringent packaging needs of medical device manufacturers. Our solutions are ideal for blister packaging applications when used in combination with compatible thermoform films, ensuring sterility, integrity, and ease of use. Suitable for EO & Gamma applications

- **Product Offerings**

Pre-Cut Sheets- Available in customized shapes with specialized round corner cuts, enabling a smooth peel and enhanced user convenience.

- **Roll Form**

Supplied in 200-meter rolls for automated form-fill-seal (FFS) operations.

- **Raw Material Options**

- **Coated HDPE** – Ensures high seal integrity with breathable properties.

Coated Hypak® – Recognized globally for superior microbial barrier and clean peel performance.



Raw Material Options

- **Easy Peel Design** – Pre-cut shapes with round corners allow clean, fiber-free opening, ideal for sterile field applications.
- **Custom Sizing Available** – Tailored dimensions and cuts to fit your packaging equipment and tray configurations.
- **Consistent Seal Performance** – Compatible with industry-standard medical-grade films for reliable sealing.
- **Regulatory Compliance** – Produced in facilities aligned with ISO 11607 and cleanroom manufacturing protocols.
- **Versatile Use** – Ideal for a wide range of medical devices including syringes, catheters, diagnostic kits, and surgical instruments.

STERILITY MONITORING



BOWIE DICK TEST PACK

Bowie & Dick test Pack for air-removal test are manufactured as per the EN ISO 11140-1:2014 compliant with EN ISO 11140-4

- Used to detect air leakages & residual condensate in Autoclave machine
- Precise and sharp indicator color transition from blue to pink change provides distinctive pass/fail result
- Each test pack has Type-1 indicator dot on top for easy identification of processed or unprocessed packages
- 100%lead free and heavy metal free indicator ink



PRODUCT CODE	PRODUCT NAME	PACKAGING
RSMBD-121/8.3	BOWIE DICK TEST PACK TYPE-2 121°C-8.3MIN/134°C-3.5MIN	120 Pcs/Ctn
RSMBD-121/10	BOWIE DICK TEST PACK TYPE-2 121°C-10MIN134°C-3.5MIN	
RSMBD-121/11	BOWIE DICK TEST PACK TYPE-2 121°C-11MIN134°C-3.5MIN	
RSMBD-121/15	BOWIE DICK TEST PACK TYPE-2 121°C-15MIN134°C-3.5MIN	
RSMBDS	BOWIE DICK TEST SHEET-A4	50 Pcs/Bag

HELIX TEST STRIPS

- Ensures air removal and steam penetration in pre-vacuum sterilizers.
- EN ISO 11140-1:2014 and EN ISO 11140-4
- Precise and sharp process indicator colour change pink to dark brown/black indicates sucessful sterilization

PRODUCT CODE	PRODUCT NAME	PACKAGING
RBDHS4S	Helix Strip-4	10 Strips/Sheet
RBDHS6S	Helix Strip-6	25 Sheets/Box



STERILITY MONITORING




CHEMICAL INDICATOR STRIPS

Chemical Indicator strips are essential tools for monitoring sterilization. They provide quick, reliable visual confirmation for the critical parameters like temperature, pressure, or chemical exposure have been met. EN ISO 11140-1:2014 compliant

- Self adhesive strips
- For Steam Application available in Type 4,5&6
- For EO application available in Type 5 which ensures accurate monitoring in EO sterilizers
- For Plasma application available in Type 4, confirms effective plasma sterilization conditions
- 100% lead free and heavy metal free indicator ink

PRODUCT NAME	PRODUCT CODE
ORO CHEMICAL INDICATOR STEAM TYPE-4 134°C-3.5MIN/121°C-15min	RSMCIST4-121/15
ORO CHEMICAL INDICATOR STEAM TYPE-4 134°C-7MIN/121°C-20min	RSMCIST4-121/20
ORO CHEMICAL INDICATOR STEAM TYPE-5 134°C-3.5MIN/121°C-15MIN	RSMCIST5-121/15
ORO CHEMICAL INDICATOR STEAM TYPE-5 134°C-7MIN/121°C-20MIN	RSMCIST5-121/20
ORO CHEMICAL INDICATOR STEAM TYPE-6 134°C- 3.5MIN/121°C-15MIN	RSMCIST6-121/20
ORO CHEMICAL INDICATOR STEAM TYPE-4 134°C-7MIN/121°C-20min	RSMCIST6-121/15
ORO CHEMICAL INDICATOR EO TYPE-5	RSMCIET5
ORO CHEMICAL INDICATOR PLASMA TYPE-4	RSMCIPT4
PACKAGING	200 Strips/Box



INDICATOR TAPES

ORO provides the accurate solution for safe sterility monitoring of sterilization wraps using these indicator tapes. These are non-hazardous tapes available for all kind steam, EO and plasma sterilization methods. The sharp indicator inks give a precise color change when meeting the prescribed sterilization method. The strong adhesive provides a firm & secure sealing keeping the integrity of the package intact, giving a comfortable & clean peeling leaving no residues behind while unpacking the package in aseptic area.

- High quality non-toxic adhesive
- Compatible to use over SMS wrap sheets, crepe sheets paper & film
- Available in 19mm & 25mm width with 50m length



CODE	RSMSIT	RSMEIT	RSPIT
TYPE	Steam	EO	Plasma
PKG	48 Pcs/Ctn		

PROCESS INDICATOR LABEL

ORO 3-line type-1 process indicator label on exposure to the sterilization process, helping to differentiate between processed and unprocessed loads

- 100% lead free and heavy metal free indicator ink
- Precise & sharp process indicator colour change
- High quality non toxic adhesive
- Packaging 1 roll x 520 nos. labels
- ISO 11140-1 compliant

CODE	RPILS	RPILE	RPILP
TYPE	Steam	EO	Plasma
PKG	60 Rolls / Ctn		



*For Label Gun please contact ORO sales representative

SHOE COVER

- The ORO universal size Shoe Cover made of Non-Woven or LDPE Fabric.
- It is water resistant, anti-static & Anti-Skid.
- The shoe Cover gives great comfort for long hours working with great comfort & extreme hygiene.



PRODUCT CODE	PRODUCT NAME	PACKAGING	
RSHCN	SHOE COVER NON-WOVEN	100 Pcs / Bag	28 Bag / Ctn
RSHCL	SHOE COVER LDPE	100 Pcs / Bag	50 Bag / Ctn

SMS JACKET

SMS Jackets are lightweight, breathable protective garments made from Spunbond-Meltblown-Spunbond material. They provide excellent barrier protection against fluids and contaminants while ensuring comfort and durability, making them ideal for healthcare use.

- Knitted cuff & caw for comfort and complete enclosure
- Soft breathable & fluid resistant fabric assuring all day protection
- Pockets & snap closure



CODE	PRODUCT NAME	PACKAGING
RJSB	SMALL, MEDIUM, LARGE, XL	10 Pcs / Bag 10 Bags / Ctn
RJSL	SMALL, MEDIUM, LARGE, XL	
RJSP	SMALL, MEDIUM, LARGE, XL	

SURGEON CAP

The disposable surgical caps are meant to be used by doctors and nurses while conducting surgeries. ORO caps are manufactured using optimum quality fabrics, which gives great comfort & highly breathable. Soft fabric with smooth texture makes the caps light weighted Backside elastic for easy size adjustment.



PRODUCT CODE	PRODUCT NAME	COLOUR	PACKAGING
RSCNW	SURGEON CAP	BLUE	100 Pcs / Bag

BOUFFANT CAP

ORO's bouffant caps provide the perfect fit to provide full hair coverage, and to stop cross-contamination by hair. These caps are highly breathable, light weight and antistatic, with an elegant look.



PRODUCT CODE	PRODUCT NAME	TYPE	PACKAGING
RBCS18B, RBCS21B	BOUFFANT CAP BLUE	STRIP 18" & 21"	100 Pcs / Bag
RBCS18W, RBCS21W	BOUFFANT CAP WHITE	STRIP 18" & 21"	100 Pcs / Bag

SAFETY WEAR



SMS LABCOAT

ORO's Lab Coats are considered as elegant, professional wear, which also give better protection. The soft, breathable and fluid-resistant material offers matchless comfort.

- Made of non-woven SMS
- Knitted cuff & collar for comfort and complete enclosure
- Soft breathable & fluid resistant fabric assuring all day protection
- Pockets & snap closure
- Available sizes- S,M,L, XL, XXL

PRODUCT CODE	SIZE	PACKAGING
RLCSB	S, M , L , XL, XXL	100 Pcs / Bag

*Customized sizes available on demand



*Please check for MOQ of colours except blue

ISOLATION GOWN, STERILE SURGEON GOWN & REINFORCE GOWN

The non-woven isolation gowns are highly breathable providing best protection

- Excellent fluid resistant nonwoven fabric
- High microbial barrier
- Great breathability giving all day comfort
- Ultimate protection against dirt & dry particulates
- Ideal for sterile area usage

CODE	TYPE	PACKAGING
RIGR	REINFORCED STERILE	10 Pcs/ Bag
RIG	ISOLATION GOWN	



FACE MASK

ORO non-absorbent & non-irritating masks remains cool on the skin providing ultimate comfort

- Fluid Resistant mask
- Adjustable nose strip
- Soft ear loop for comfortable wear
- TYPE II IR mask conforming to EN14683:2019+AC:2019
- $\geq 98.4\%$ BFE, Superior Melt blown
- Ear loop & Tie-on variants available
- Available in 3PLY & 4PLY



*Colour available for 3 Ply & 4 Ply



CODE	PRODUCT NAME	PACKAGING
RFM3E	FACE MASK 3 PLY EAR LOOP (BLUE, BLACK, PINK, WHITE)	50 Pcs/ box
RFM4E	FACE MASK 4 PLY EAR LOOP (BLUE, BLACK PINK, WHITE)	50 Pcs/ box
RFM3T	FACE MASK 3 PLY TIE ON (BLUE, BLACK PINK, WHITE)	50 Pcs / box
RFM4T	FACE MASK 3 PLY TIE ON (BLUE, BLACK PINK, WHITE)	50 Pcs / box

DISPOSABLE COVERALL

Disposable Coverall is designed to cover the entire body and avoid cross-contaminations in the cleanroom environment and process carried out in the cleanroom. Single use disposable coverall act as a barrier for cleanroom contamination by human activities.

- Suitable for controlled environment
- Pre-sterile for quick and easy access
- Made of breathable PP+PE Non-Woven Fabric
- EN13795-2 complaint

*For details contact ORO sales representative



REUSABLE COVERALL

Reusable Coverall plays key role in the contamination control of clean room environment. Human is the sole source of contamination in clean room environments apart from the material. Human skin carries various types of microorganisms like Bacteria, Fungal spores and viruses which can be carried with the scale flax and contaminate the clean room environment. ORO Reusable Coverall are manufactured by using 98% Polyester and 2% carbon Fibers, which have antistatic properties. Reusable Coverall are designed and packed in such a way that aseptic operations can be fulfilled. 3 Piece Reusable coverall covers complete body from top to toe and avoid the contaminations through particles shading in the cleanroom. ORO Reusable Coverall are breathable to achieve comfortable operations.

*For details contact ORO sales representative





Reach Global India Pvt. Ltd.

- **Reach Global India Private Limited.**
Office No. 1&2, Survey No. 64/1, Wing-A, First Floor,
Sai Raj Regency, Pimple Saudagar, Pune: 411027
- **Plot No 1-22, Gat No. 1064, Chordia Industrial Park,**
Naigaon, Shirwal, Taluka-Khandala,
Dist. Satara, State-Maharashtra (India)
Pin Code - 412801

Customer Care: +91 9370 788 850/ +91-(020) 27400987
info@orohealthcare.com

