

DENTAL CONSUMABLES & STERILIZATION PACKAGING



EN ISO 13485 | CE | MADE IN INDIA



VISION

While retaining our current reputation of "your infection guard" Our vision is to become 'leading manufacturer & supplier' of barrier packaging & associated products for controlled area.

MISSION

Our Sterilization packaging products, Protective clothing and supplementary products are targeted to cater the needs of 'Healthcare segment' for their CSSD, Production, Quality Control & Assurance sections. We aim to support their needs by offering suitable products

OUR VALUES

BELIEF

Our strong belief is our people and we work by taking ownership to deliver our promises persistently

FUTURISTIC

We plan ourselves by targeting long term and sustainable success

UNIFY

We encourage and challenge each other with respect to achieve our common objective with fun and by helping each other unified

UNDERSTANDING

We are open and have courage to try new & learn from our mistakes. we get happiness in developing partnership with our customer



REACH GLOBAL

INFECTION BARRIERS

Barrier Film.....	04
Barrier Envelope.....	04
Sleeves.....	05

SAFE STERILIZATIONS

Bowie dick Kit.....	05
Sterilization Reels.....	06
Safe Seal Pouches.....	06
Indicator Tapes.....	07
Steam Indicator Strip.....	07

IMPRESSION ACCESSORIES

Bite Tray.....	08
Fluoride Tray.....	08
Impression Tray.....	09
Mixing Bowl.....	09
Spatula.....	10
Dispensing Gun.....	10

GERM-FREE EVACUATIONS

Intra Oral Tip.....	11
Saliva Ejector.....	11
Surgical Aspirator Tip.....	12
Vented Oral Evacuation Tip.....	12

HYGIENIC DISPOSABLES

Micro Applicator, Tip.....	13
Micro Applicator, Brush.....	13
Disposable Air / Water	
Syringe Tips	14

SECURE DISPENSERS / ORGANIZERS

CompuLe Gun, Micro	14
Instrument Tray.....	15
Endo Block.....	15
Endo Boxes.....	16

RETRACTORS AND EXPANDERS

Cheek Retractor.....	16
Mouth Prob.....	17

PROTECTIVES & ABSORPTIVES

Bouffant Cap.....	17
Face Mask.....	18
Face Shield.....	18
Surgical Kit Dental	19
Patient Bib.....	19
Dental Apron.....	20
Protective Glass.....	20
SMS Lab Coat & Jacket.....	21
Surgical Gown.....	21
Cotton Roll & Dispenser.....	22
Sponges.....	22

SAFE KEEPERS

Bur holder.....	23
Denture Box.....	23
Retainer Box.....	24
Aligner Box.....	24

LITTLE CROWN.....27

IMPRESSION ACCESSORIES

Mixing Tip.....	28
-----------------	----

LOUPE.....29

HANDPIECE.....30-31

SURE ENDO.....32-35

INFECTION BARRIERS

BARRIER FILM



These perforated sheets with non-sticky edge are effective to reduce cross contamination with single-use and covers anything you might touch during procedure. Easy to protect the areas which are difficult to clean.

Smart packaging with self-dispense is best for direct usage with easy tear off.

SPECIFICATIONS

Code	Size	Packaging
Blue, RBFB	4" X 6", 1200 Sheets	With stand-1 wrap / box, 8 boxes / case
Clear, RBFC	4" X 6", 1200 Sheets	Without stand-1 wrap / box, 12 boxes / case

INFECTION BARRIERS

DISPOSABLE BARRIER ENVELOPES



Barrier Envelope protects the phosphor storage plates from mechanical damage and information loss caused by ambient light. The Envelope keeps the phosphor storage plate clean while reducing cross contamination.

Peel and stick end provide a consistent seal and tear off.

SPECIFICATIONS

Code	Size	Packaging
EBE0	Size 0	100 Pcs / box, 100 boxes / case
RBE1	Size 1	100 Pcs / box, 100 boxes / case
RBE2	Size 2	100 Pcs / box, 100 boxes / case

INFECTION BARRIERS

SLEEVES



Sleeves are important to maintain hygienic workflow and are vital to operatory asepsis, which play crucial role in infection control by reducing clean up time and avoid cross-contamination.

We are specialized to provide sleeves for all varieties and sizes of your requirement.

SPECIFICATIONS

Type	Code	Size	Packaging
Sensor Sleeve	RSLS8X16	8X16	500 pcs/box, 36 boxes / case
Scaler Sleeve	RSLSC1.4X7	1.4X7	500 pcs/box, 36 boxes / case
Intra Oral Camera Sleeve	RSLIC10.2X2	10.2X2	500 pcs/box, 40 boxes / case
Light Cure Sleeve	RSLLC4.3X2.95	4.3X2.95	200 pcs/box, 42 boxes / case
OPG Bite Block Sleeve	RSLOBB2.7X1.1	2.7"X1.1"	500 Pcs/Box,
Light cure Sleeve Large	RSLLC10.2x2.2	10.2" x2.2"	500 pcs/box, 40 boxes / case

SAFE STERILIZATIONS

BOWIE DICK



The test pack consists of a series of air removal and steam penetration barriers. Single use test packs are designed for validation in sterilization process in the sterilization process.

- Colour change from Blue to Black will be clear and stable, and will not disappear easily
- Manufactured in the absence of lead or other heavy metals
- Design to spot any malfunction / lack of efficacy in the pre-vacuum sterilizer

SPECIFICATIONS

Code	Type	Pkg
RSMBD-121/8.3	121°C-8.3min/134°C-3.5MIN	120 Pcs/ Ctn
RSMBD-121/10	121°C-10min/134°C-3.5MIN	
RSMBD-121/11	121°C-11min/134°C-3.5MIN	
RSMBD-121/15	121°C-15min/134°C-3.5MIN	
RSMBD	ORO BD test sheet-A4	

SAFE STERILIZATIONS

STERILIZATION REELS

CE CERTIFIED



Tripled sealed flat rolls are constructed from transparent PET/CPP multi-layers copolymer film and medical grade paper. The process indicators for steam and ethylene oxide are water-based, non-toxic inks which provide accurate colour change indication to confirm the process is achieved.

SPECIFICATIONS

Code	Size	Packaging
RRM50	50mm x 200mm	10 rolls / case
RRM55	55mm x 200mm	10 rolls / case
RRM75	75mm x 200mm	8 rolls / case
RRM100	100mm x 200mm	6 rolls / case
RRM125	125mm x 200mm	8 rolls / case
RRM150	150mm x 200mm	4 rolls / case
RRM200	200mm x 200mm	4 rolls / case
RRM250	250mm x 200mm	4 rolls / case
RRM300	300mm x 200mm	2 rolls / case
RRM350	350mm x 200mm	2 rolls / case
RRM400	400mm x 200mm	2 rolls / case
RRM450	450mm x 200mm	2 rolls / case
RRM500	500mm x 200mm	2 rolls / case

SAFE STERILIZATIONS

SELF SEAL POUCHES

CE CERTIFIED



These well-designed sterilization pouches come with a triple-layered seal, delivering three levels of safety for easy, effective and infection-free sterilization. These pouches conform to the three parameters of time, temperature & pressure steam / EO valid up to 6 months. Easy monitoring can be achieved through internal and external colour change indicators.

SPECIFICATIONS

Code	Size	Packaging
RPMS 57 X 100	57mm x 100mm	200 pcs/box, 50 boxes / case
RPMS 57 X 125	57mm x 125mm	200 pcs/box, 30 boxes / case
RPMS 90 X 160	90mm x 160mm	200 pcs/box, 30 boxes / case
RPMS 70 X 255	70mm x 255mm	200 pcs/box, 20 boxes / case
RPMS 90 X 260	90mm x 260mm	200 pcs/box, 20 boxes / case
RPMS 100 X 280	100mm x 280mm	200 pcs/box, 20 boxes / case
RPMS 134 X 283	134mm x 283mm	200 pcs/box, 20 boxes / case
RPMS 250 X 400	250mm x 400mm	200 pcs/box, 5 boxes / case
RPMS 300 X 400	300mm x 400mm	200 pcs/box, 5 boxes / case
RPMS 150 X 280	150mm x 280mm	200 pcs/box, 10 boxes / case
RPMS 330 X 450	330mm x 450mm	200 pcs/box, 5 boxes / case
RPMS 250 X 550	250mm x 550mm	200 pcs/box, 5 boxes / case
RPMS 190 X 360	190mm x 360mm	200 pcs/box, 10 boxes / case

*Customized sizes available on demand

SAFE STERILIZATIONS

INDICATOR TAPES



Autoclavable adhesive tapes are used to confirm whether a specific temperature has been attained. These are water-based and non-toxic chemical indicator tapes with stripes. Extra adhesive is provided to fix on all types of paper & plastic.

SPECIFICATIONS

Type	Code	Size	Packaging
STEAM	RSMSIT	19mm x 50m	48 Pcs/case
EO	RSMEIT	19mm x 50m	48 Pcs/case
PLASMA	RSMPIT	19mm x 50m	48 Pcs/case

SAFE STERILIZATIONS

STEAM INDICATOR STRIPS



Steam Indicator Strips are designed to be used in steam sterilizers operating at either 121°C for 15 minutes or 134°C for 3-5 minutes. The water-based non-toxic chemical indicator will show a clear and accurate colour change from blue to black after successful steam sterilization.

SPECIFICATIONS

Name	Type	Code	Packaging
Steam Indicator Strip	4	RSMCIST4-121/15	200 Strips/Box
		RSMCIST4-121/20	200 Strips/Box
	5	RSMCIST5-121/15	200 Strips/Box
		RSMCIST5-121/20	200 Strips/Box
	6	RSMCIST6-121/15	200 Strips/Box
		RSMCIST6-121/20	200 Strips/Box
EO Indicator Strip	4	RSMCIET4	200 Strips/Box
	5	RSMCIET5	200 Strips/Box
Plasma Indicator Strip	5	RSMCIPT4	200 Strips/Box

IMPRESSION ACCESSORIES

BITE TRAY



Bite registration tray is an easy and economical way to reduce the risk of cross-contamination and infection. The robust construction of the single use dual-bite impression tray is made with polished plastic, smooth edges to ensure the patient's comfort. The net gives best holding capability to under-lock the material, and retains it in the same position.

SPECIFICATIONS

Code	Type	Packaging
RBTAT	Anterior	35 pcs/bag, 20 bag/case
RBTAS	Assorted	50 pcs/bag, 20 bag/case
RBTFA	Full Arch	35 pcs/bag, 20 bag/case
RBTPT	Posterior	50 pcs/bag, 20 bag/case
RBTQD	Quadrant	30 pcs/bag, 20 bag/case
RBTSL	Sideless	35 pcs/bag, 20 bag/case

IMPRESSION ACCESSORIES

FLUORIDE TRAY




Disposable fluoride trays available in a variety of sizes are custom-fitted foam trays for fluoride gel, to control decay and strengthen teeth structure. They are designed to offer a natural, contoured fit and exceptional comfort for fluoride treatment.

SPECIFICATIONS

Code	Size & Colour	Packaging
RFTL	Large, Blue	50 pcs/bag, 50 bags / case
RFTM	Medium, Yellow	50 pcs/bag, 60 bags / case
RFTS	Small, Green	50 pcs/bag, 80 bags/case

IMPRESSION ACCESSORIES

IMPRESSION TRAY



Autoclavable

Disposable

Perforated impression trays are made from strong, rigid materials to generate consistent impressions. These trays are distortion-free, pre-shaped, readily trimmed, and are flame-softened for adjustment. The handle features non-slip grip and are compatible with wide variety of impression materials. Available in disposable and autoclavable material up to 121°C.

SPECIFICATIONS

Code	Size	Packaging
RIMTA1&2	Size 1&2	10 pcs/bag, 50 bags / case
RIMTA3&4	Size 3&4	10 pcs/bag, 50 bags / case
RIMTA5&6	Size 5&6	10 pcs/bag, 50 bags / case
RIMTA7&8	Size 7&8	10 pcs/bag, 50 bags / case
RIMTA9&10	Size 9&10	10 pcs/bag, 50 bags / case
RIMTA1-10	Size 1-10	10 pcs/bag, 50 bags / case

IMPRESSION ACCESSORIES

MIXING BOWL



Flexible PVC bowl work that is best for mixing impression material and is easy to re-use after cleaning. The material is suitable for both hard and soft mixtures.

SPECIFICATIONS

Code	Capacity	Colour	Packaging
RMBL600	600ml	Blue	100 pcs/case
RMBM500	500ml	Green	100 pcs/case
RMBS280	280ml	Clear Green/Red	180 pcs/case
RMBXS150	150ml	Clear Green/Red	250 pcs/case

IMPRESSION ACCESSORIES

SPATULA



A flat and smooth fore-end spatula is good for mixing. Autoclavable up to 121°C.

SPECIFICATIONS

Code	Colour	Packaging
RSP	Pink, White, Yellow & Blue	50 pcs/bag, 10 bags / case

IMPRESSION ACCESSORIES

DISPENSING GUN



Dispensing gun works perfectly to dispense dental adhesive silicones and dental materials for perfect registration of dental impression. This sturdy device is made of high-quality plastic and consists of a release lever with a plunger, which coordinates simultaneously to load and remove the cartridges. Works perfectly with ORO mixing tips.

SPECIFICATIONS

Code	Type	Packaging
RDG1:1X50ML	1:1/1:2 50 ML.	1 pc/box, 50 boxes / case
RDG1:1X25ML	1:1/1:4 25 ML	1 pc/box, 50 boxes / case
RDG1:10X50ML	1:10/1:4 50 ML	1 pc/box, 50 boxes / case

GERM FREE EVACUATIONS

INTRA ORAL TIP



Intra oral tip description : Oro Intra Oral Tips are dispensing tips, especially designed for controlled release of core buildup material and other medium viscosity flowable materials.

SPECIFICATIONS

Code	Type	Packaging
RIOTN1	N1	100 pcs/bag
RIOTN2	N2	100 pcs/bag
RIOTN3	N3	100 pcs/bag
RIOTN5	N5	100 pcs/bag
RIOTN8	N8	100 pcs/bag

GERM FREE EVACUATIONS

SALIVA EJECTOR



The non-sterile disposable saliva ejector is made of non-toxic material. The incorporated copper wire support to maintain a stable and precise shape. The non-removable and special shape tip guarantees optimal suction, preventing any damage to the Patient's oral mucosa.

SPECIFICATIONS

Code	Type	Packaging
RSEB	Blue	100 pcs/bag, 50 bags/case
RSEC	Clear	100 pcs/bag, 50 bags/case
RSEG	Green	100 pcs/bag, 50 bags / case

SURGICAL ASPIRATOR TIP



ORO's surgical aspirator tips are made of high-quality polypropylene and designed to fit all high-volume evacuations. These tips are rigid enough to aspirate and retract the cheek without bending, and the 30° permanent angle allows better visual contact with the site, while the smooth edges enhance patient comfort. These can be disposed after use or autoclaved at 121°C.

SPECIFICATIONS

Code	Type	Packaging
RSPT1/16B	Blue	100 pcs/bag, 50 bags/case
RSPT1/14G	Clear	100 pcs/bag, 50 bags/case
RSPT1/18W	Green	100 pcs/bag, 50 bags / case

VENTED ORAL EVACUATION TIP



ORO's single use disposable tips contain smooth edges for maximum patient comfort. One end is vented and the other is non-vented. These tips fit all standard evacuator hose ends and are non-sterilized.

SPECIFICATIONS

Code	Size & Colour	Packaging
RVOET6	6", white	100 pcs/bag, 10 bags/case
RVOET3	3", white	100 pcs/bag, 20 bags / case

HYGIENIC AND DISPOSABLES

MICRO APPLICATOR TIP

CE CERTIFIED



Micro applicator tip is non-linting and non-absorbent with a featured bendable handle. These tips are designed for precise placement of bonding agents, etchants, haemostatic solutions, cavity liners and disclosing solutions. They facilitate the proper placement of materials into difficult areas. Available in multiple colours for easy identification in a variety of sizes.

SPECIFICATIONS

Code	Colour	Packaging
RMAFW	Superline White	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMAFM	Fine Mauve	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMAFY	Fine Yellow	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMARB	Regular Blue	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMARG	Regular Green	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMARO	Regular Orange	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMARP	Regular Purple	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMA SBL	X-Slim, Black	100 pcs/tube, 4 tubes/box, 32 boxes/case

HYGIENIC AND DISPOSABLES

MICRO APPLICATOR BRUSH




Micro applicator brush is non-linting and non-absorbent with a featured bendable handle. These brushes are designed for precise placement of bonding agents, etchants, haemostatic solutions, cavity liners and disclosing solutions. They facilitate the proper placement of materials into difficult areas. Available in multiple colours for easy identification in a variety of sizes.

SPECIFICATIONS

Code	Colour	Packaging
RMABW	White	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMABG	Green	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMABBL	Blue	100 pcs/tube, 4 tubes/box, 32 boxes/case
RMABY	Yellow	100 pcs/tube, 4 tubes/box, 32 boxes/case

DISPOSABLE AIR/WATER SYRINGE TIPS



These tips have an inner core water passage channel that fits nicely into the syringe preventing water contamination in the air duct. The core design of the tips provides instant dry air and prevents cross-contamination while delivering superior air / water spray ratios.

SPECIFICATIONS

Code	Colour	Packaging
RAWSDA	Assorted	250 pcs / box, 50 boxes / case
RAWSDB	Blue	250 pcs / box, 50 boxes / case
RAWSDC	Clear	250 pcs / box, 50 boxes / case

COMPULE GUN



ORO's compule gun is designed to deliver the material with easy handling of the compule tips, i.e. direct and precise placement of the composite and compomer restorative into the cavity. These injection guns match most charging barrels, and are the best choice for sticky materials.

SPECIFICATIONS

Code	Colour	Packaging
RCMG	150mm x 100mm	1 pcs / bag, 160 pcs / case

SECURE DISPENSER / ORGANIZER

INSTRUMENT TRAY

Disposable



Autoclavable



Stylish cup dispensers are a smart way to keep your clinic and office hygienic and organized. ORO provides a variety of dispensers like wall mounted or chair mounted for different purposes. The cups are sturdy and carry 5 OZ / 148ml, and smartly fit into ORO dispensers. Out tissue box bring handiness for easy single tissue pull-out, while keeping the remaining tissues safe.

SPECIFICATIONS

Code	Type	Packaging
RITA22X16.5, RITA33X24, RITA20X10	Autoclavable	160 pcs/case
RITD28X19B-Blue	Disposable	100 pcs/bag 4 bag / case
RITD28X19G-Green		
RITD28X19O-Orange		
RITD28X19W-White		
RITD28X19Y-Yellow		

SECURE DISPENSER / ORGANIZER

ENDO BLOCK



Graduated holes of different lengths provide an accurate, quick method of adjusting root canal instruments at the required length.

SPECIFICATIONS

Code	Colour	Packaging
REB	Blue, Grey	250 pcs / box

ENDO BOXES



These organizers are made of plastic and have removable, colour-coded silicon cups for Gutta Percha points & Paper point along with storage space for burs / files / reamers / spreaders. Autoclavable up to 135°C

SPECIFICATIONS

Code	Shape	Packaging
REB36	Square, 36 holes	1 PC / box, 120 boxes / case
REB30	Circular, 30 holes	1 PC / box, 64 boxes / case
REB60	Rectangle, 60 holes	1 PC / box, 96 boxes / case

CHEEK RETRACTOR



This accessory is for retraction of lips and cheek during a dental treatment. Cheek Retractors are durable, flexible and convenient to use. The smooth surface makes the patient comfortable.

SPECIFICATIONS

Code	Type	Packaging	Autoclave
RCRCL, RCRCLM, RCRCS	C type: Large, Medium, Small	2 pcs / box, 200 boxes / case	Upto 134°C
RCRML, RCRMM, RCRMS	M type: Large, Medium, Small	2 pcs / box, 200 boxes / case	Upto 134°C
RCRYL, RCRLS	Y type: Large, Medium, Small	2 pcs / box, 200 boxes / case	Not Applicable

RETRACTORS AND EXPANDERS

MOUTH PROP



The upper and lower teeth can be easily propped up in the mouth. These props are made from comfortable and elastic medical plastic material, with no glue and peculiar smell. These props are helpful to avoid patients' fatigue. The tongue guards block tongue movement to eliminate the chance of accidental injury. These props are resistance to high temperatures of up to 121°C, and can also withstand high-pressure disinfection.

SPECIFICATIONS

Code	Size	Packaging
RMPL	Large	50 pcs / bag, 500 bags / case
RMPM	Medium	50 pcs / bag, 500 bags / case
RMPS	Small	50 pcs / bag, 500 bags / case

PROTECTIVES & ABSORPTIVE

BOUFFANT CAP



ORO's bouffant caps provide the perfect fit to provide full hair coverage, and to stop cross-contamination by hair. These caps are highly breathable, light weight and antistatic, with an elegant look.

SPECIFICATIONS

Code	Type	Packaging
RBCS18B	Strip 18" Blue	100 pcs / bag, 50 bags / case
RBCS18B	Strip 18" White	100 pcs / bag, 50 bags / case
RBCS21B	Strip 21" Blue	100 pcs / bag, 50 bags / case
RBCS21W	Strip 21" White	100 pcs / bag, 50 bags / case
RSCNW	BLUE (SURGEON CAP)	100 pcs / bag, 50 bags / case

FACE MASK



ORO's 3-ply gentle, non-absorbent and non-irritating mask remains cool on the skin for ultimate ease and comfort. The adjustable nose piece helps to customize the fit, while the ear loop ties are soft, stretchable and latex free. Follow size and BFE as 17.5cm X 9.5cm, BFE $\leq$ 95%–99.05%.

SPECIFICATIONS

Code	Type	Packaging
RFM3EB	3 Ply Ear Loop Blue	50 pcs / box, 50 boxes / case
RFM3EP	3 Ply Ear Loop Pink	50 pcs / box, 50 boxes / case
RFM3EW	3 Ply Ear Loop White	50 pcs / box, 50 boxes / case
RFM3EBL	3 Ply Ear Loop Black	50 pcs / box, 50 boxes / case
RFM3TB	3 Ply Tie on Blue	50 pcs / box, 40 boxes / case
RFM4EBL	4 Ply Ear Loop Black	50 pcs / box, 50 boxes / case
RFM4EP	4 Ply Ear Loop Pink	50 pcs / box, 50 boxes / case
RFM4EW	4 Ply Ear Loop White	50 pcs / box, 50 boxes / case
RFM4EB	4 Ply Ear Loop Blue	50 pcs / box, 50 boxes / case

FACE SHIELD



Anti-fog full face protection shield protects the face and eyes from flying debris and blood / water splashes.

- The foam shield is disposable and can be used over eyeglasses. The fit is very comfortable with elastic.
- The re-usable shield is washable, easy to clean, and light weight for better comfort. The shields take quick position and fits over eyeglasses.

SPECIFICATIONS

Code	Type	Packaging
RFSD33*32	Face Shield Disposable 33x32 cm	10 pcs/bag, 20 bags / case
RF5FX10S	With frame & 10 shields	1 set / box, 30 boxes / case

PROTECTIVES & ABSORPTIVE

SAFETY WEAR SURGICAL KIT DENTAL



Implant kits are designed to maximize efficiency and simplify surgical procedures. We offer a range of surgical kits specially configured to meet doctors' needs for virtually all clinical scenarios. The kit layout features easy-to-follow protocol flow, to assist in economizing time and effort for every phase of surgery.

SPECIFICATIONS

Code	Size	Packaging
RSKD	Instrument table cover 80cmx100cm	20 Kits/ case
	3 ply Earloop face mask	
	Bouffant Cap	
	Tie back surgical gown	
	Shoe Cover	
	Towel 25 cmx25cm	
	U non-woven drape with adhesive (140cmx170cm)	
	Non-woven sponges 3x3x4ply	
	Tubing sleeve (plastic sheet) with adhesive 8.5cmx120cm)	

PROTECTIVES & ABSORPTIVE

PATIENT BIB



ORO's 3-ply patient bib is made of soft tissue and poly of highly absorbent nature. The bonding of 2-ply tissues and 1-ply poly work securely together to eliminate separation of layers, and the reinforced edges provide added strength and durability. The tissues resist ripping even when wet, and can therefore double as a counter top drape for soaking up fluids.

Twisted clips are suitable to be used with the bibs while operating on the patient. The sufficient length ensures comfort.

SPECIFICATIONS

Code	Colour	Packaging
RPBB	Blue	125 pcs / bag, 4 bags / case
RPBG	Green	125 pcs / bag, 4 bags / case
RPBL	Lavender	125 pcs / bag, 4 bags / case
RPBO	Orange	125 pcs / bag, 4 bags / case
RPBP	Pink	125 pcs / bag, 4 bags / case
RPBW	White	125 pcs / bag, 4 bags / case
RPBY	Yellow	125 pcs / bag, 4 bags / case

PROTECTIVES & ABSORPTIVE

DENTAL APRON



Dental apron help to protect skin and clothing from damage such as liquid splashing in the torso area. The apron has good absorbing capacity and its pattern is perfect to channelize liquid downward.

SPECIFICATIONS

Code	Type	Packaging
RPANWTB	Non woven Tie Back	50 pcs / bag, 25 bags /case
RPAVTB	Vinyl Tie Back	50 pcs / bag, 25 bags /case
RPANWRN	Round Neck	50 pcs / bag, 25 bags /case

PROTECTIVES & ABSORPTIVE

PROTECTIVE GLASS



Protective glasses are specially designed to protect the eyes from injury. These glasses are anti-fog, anti-hard light, and anti-ultraviolet ray. With a smooth surface they provide great comfort. The glasses also feature an anti-scratch material, which offers excellent protection and opal view.

SPECIFICATIONS

Code	Type	Packaging
REWO	Orange	12 pcs / box, 12 boxes / case
REWC	Clear	12 pcs / box, 12 boxes / case

PROTECTIVES & ABSORPTIVE

SMS LAB COAT & JACKET



SMS protective lab coats and jackets come with pockets and knitted wrist cuffs. These are considered as elegant, professional wear, which also give better protection. The soft, breathable and fluid-resistant material offers matchless comfort.

SPECIFICATIONS

Product	Code	Size	Colour	Pkg
Labcoat	RLCSBXXL, RLCSBXL, RLCSBL, RLCSBM, RLCSBS	XXL, XL, L, M, S	Blue, Pink	10 pcs / bag, 10 bags / case
Jacket	RJSBL, RJSBM, RJSBS	L, M, S	Blue, Pink	10 pcs / bag, 10 bags / case

PROTECTIVES & ABSORPTIVE

ISOLATION & SURGEON GOWN



ORO's economical, disposable barrier gowns are made of fluid-repelling spun bond polypropylene. These gowns offer basic protection against light fluid spray, dirt, and dry particulates. This popular fabric is economical, breathable, light weight and flexible. The multi-purpose isolation gown features elasticized wrists, neck and waist ties, and overlapping back panels.

SPECIFICATIONS

Code	Description	Packaging
RIGNS	Isolation Gown non-sterile	10 pcs / bag, 10 bags / case
RIGS	Surgeon gown sterile	1 pc / bag, 50 bags / case

PROTECTIVES & ABSORPTIVE

COTTON ROLL AND DISPENSER



Cotton rolls are highly absorbent and made using 100% pure virgin white cotton with long fibres. These rolls are highly flexible, conform easily and maintain their shape during use. The absorption capacity is 10 times of its normal weight.

A single wrap cotton roll fits perfectly in the dispenser which gives it a smart look, while maintaining safety from cross-contamination. It has a rubber base for a sturdy non-slip surface grip with a transparent body.

SPECIFICATIONS

Product	Code	Size	Packaging
Cotton Roll	RCR1X500	#1, 1-1/2 in X 5/16 in	500 pcs / bag 80 bags / case
Cotton Roll	RCR2X500	#2, 1-1/2 in X 3/8 in	500 pcs / bag 70 bags / case
Cotton Roll	RCR3X500	#3, 1-1/2 in X 1/2 in	500 pcs / bag 60 bags / case
Cotton Roll Dispenser	RCDP	120mm x 55mm x 97mm	1 pcs / bag 40 bags / case
Cotton Roll	RCR1X1000	#1, 1-2/2 in X 5/16 in	1000 pcs / bag 40 bags / case
Cotton Roll	RCR2X1000	#2, 1-1/2 in X 3/18 in	1000 pcs / bag 35 bags / case
Cotton Roll	RCR3X1000	#3, 1-1/2 in X 1/2 in	1000 pcs / bag 30 bags / case

PROTECTIVES & ABSORPTIVE

SPONGES



ORO's non-sterilized non-woven sponges are soft, smooth and non-linting. The material offer maximum absorbency and little adhesion to wounds with pure white quality. Soft and comfortable to use either dry or moistened, these sponges absorb 10 times their weight.

SPECIFICATIONS

Code	Size	Packaging
RSPGN4P2X2	4Ply, 2"x2"	200 pcs / Sleeve, 60 Sleeves / case
RSPGN4P3X3	4Ply, 3"x3"	200 pcs / Sleeve, 24 Sleeves / case
RSPGN4P4X4	4ply, 4"x4"	200 pcs / Sleeve, 24 Sleeves / case

SAFE KEEPERS

BUR HOLDER



Bur holders will organize your instruments for use during the procedure and while sterilization process. You can use these holders during procedures to give you quick and easy accessibility to your burs. Varying hole sizes can accommodate different types of shanks and cutting heads.

SPECIFICATIONS

Product	Shape & Shot	Packaging
RHBHAC60	Circular, 60 holes, Autoclavable	1 pc / box, 60 boxes / case
RBHAS24	Square, 24 holes, Autoclavable	1 pc / box, 100 boxes / case

SAFE KEEPERS

DENTURE BOX



ORO's denture box is a robust storage container for hygienic storage of dentures, braces, gum shields and other orthodontic appliances. The boxes are made of superior plastic to keep them healthy and clean.

SPECIFICATIONS

Code	Description	Colour	Packaging
RDBB RDBG RDBLB RDBR	Denture Box	Blue & Green, Light Blue, Red	12 pcs/box 16 boxes / case 296 pcs /case
RDBT	Denture Box with Tray	Dark Green, Green with white box	180 pcs / case

SAFE KEEPERS

RETAINER BOX / ALIGNER BOX



ORO's retainer box is made of sturdy plastic to protect the patient's retainer from bacteria or breakage. The retainer box includes vent holes or a hinged lid.

SPECIFICATIONS

Code	Colour
Square RRBS	Yellow, Clear, Blue, Clear Green, Rose, White red,
Triangle RRBT	Yellow, Blue, Green, Rose Red, White
RABY	Yellow
RABG	Green
RABP	Pink
RABPL	Purple
RABB	Blue
RABW	White
RABB	Black

MISCELLANEOUS

GLASS DAPPEN DISH



These are dense glass vessels for mixing dental medicaments or fillings. The thick and heavy base makes the dish stable and offers protection from spill. It is ideal for mixing acrylics, amalgam, resins, sealants and other materials.

SPECIFICATIONS

Code	Colour	Packaging
RDDG	Amber, Blue, Green, Clear	12 pcs / box, 18 boxes / case

MISCELLANEOUS

PRE-BENT TIPS



Pre-bent tips are disposable, blunt-end needle tips designed to dispense etchants, resins and flowable composites cleanly and easily. These tips are offered in several different gauge sizes to accommodate most materials and procedures.

SPECIFICATIONS

Code	Gauge	Pkg
RPBT18	18ga (cavity liners, bases and composites)	100 pcs / bag
RPBT19	19ga (flowable composites and sealants)	100 pcs / bag
RPBT20	20ga (flowable composites)	100 pcs / bag
RPBT22	22ga (sealants and cements)	100 pcs / bag
RPBT25	25ga (etchants)	100 pcs / bag

MISCELLANEOUS

WOODEN WEDGES



Eco-friendly wooden wedges are anatomically shaped for better and safer procedures. We offer 6 wedges which can be trimmed without splintering, for a perfect custom fit. The wedges' colour-coding agent does not recolour* the tooth or restoration upon removal.

SPECIFICATIONS

Code	Colour	Packaging
RWWA	Assorted	100 pcs / box, 36 boxes / case
RWWB	Blue	100 pcs / box, 36 boxes / case
RWWG	Green	100 pcs / box, 36 boxes / case
RWWO	Orange	100 pcs / box, 36 boxes / case
RWWP	Pink	100 pcs / box, 36 boxes / case
RWWW	White	100 pcs / box, 36 boxes / case
RWWY	Yellow	100 pcs / box, 36 boxes / case

PROPHY BRUSH & CUP



ORO's prophy brush reaches into recesses where a normal brush can't reach. The variant shape of the brush removes heavy stains from occlusal surfaces and hard-to-reach inter-proximal spaces. Prepares surface for sealant and orthodontic procedure. It is made up of nylon, and is latex-free.

X-RAY FILM HOLDER

X-ray holder takes quick and precise anterior and posterior bisecting angle radiographs, and eliminates manual retention on the film packet by patient. It enhances patient comfort with smooth edges.

SPECIFICATIONS

Code	Gauge	Packaging
RXFH	Beige	4 Pcs / box, 200 boxes / case

SPECIFICATIONS

Code	Type	Image	Packaging
RPB-330	latch type flat white		100 pcs / box, 100 boxes / case
RPB-340	latch type asperred white		100 pcs / box, 100 boxes / case
RPB-350	latch type flat black		100 pcs / box, 100 boxes / case
RPB-360	latch type tapered black		100 pcs / box, 100 boxes / case
RPB-370	latch type round white		100 pcs / box, 100 boxes / case
RPB-330A	latch type flat assorted		100 pcs / box, 100 boxes / case
RPB-370A	latch type round assorted		100 pcs / box, 100 boxes / case
RPC	cup latch type style 6 rib		100 pcs / box, 100 boxes / case

*Other variants available on request.



LITTLE CROWN

Stainless Steel Primary
Molar Crown

Available in
Kit: 48/96 PCS
Refill: 5 PCS

ACCURATE OCCLUSAL ANATOMY

GINGIVAL CONTOURING

CORROSION RESISTANCE

● Crown made of Austenite alloy


● Ergonomic Design

● Accurate occlusal anatomy

● Good polishing finish

● Pre-trimmed and crimped

● Save your chair time



O R O
L I T T L E
C R O W N

*Perma Crown & Zirconia Crown are Coming Soon!

MIXING TIP



Seil Global Mixing Tips are designed for optimal efficiency in impression material application. This pack includes 50 disposable mixing ensures precise mixing and accurate application of light body impression materials, enhancing both performance and hygiene. The mixing tips are engineered to minimize waste and prevent cross contamination. Ideal for maintaining high standards of infection control and streamlining dental procedures, these tips contribute to improved workflow and material management.

MIXING TIP

SPECIFICATIONS

Code	Type	Packaging
RMTBFC1:1	Brown Flat Core 1:1	50 pcs/Bag, 50Bags/ Case
RMTBTRS1:1	Brown Tapered Resin Short 1:1	50 pcs/Bag, 50Bags/ Case
RMTBTRL1:1	Brown Tapered Resin Long 1:1	50 pcs/Bag, 50Bags/ Case
RMTYLBL1:1	Yellow Light Body Long 1:1	50 pcs/Bag, 50Bags/ Case
RMTYLS1:1	Yellow Light Body Short 1:1	50 pcs/Bag, 50Bags/ Case
RMTGHL1:1	Green Heavy Body Long 1:1	50 pcs/Bag, 50Bags/ Case
RMTGHBS1:1	Green Heavy Body Short 1:1	50 pcs/Bag, 50Bags/ Case
RMTPRVHB1:1	Pink Regular/Heavy Body 1:1	50 pcs/Bag, 50Bags/ Case
RMTTC&B1:10	Temporary Crown & Bridge 1:10	50 pcs/Bag, 50Bags/ Case
RMTBTC1:1	Blue Tapered Cement 1:1	50 pcs/Bag, 50Bags/ Case
RMTBTCL1:10	Brown Temporary Cement Long 1:10	50 pcs/Bag, 50Bags/ Case
RMTBTCS1:10	Brown Temporary Cement Short 1:10	50 pcs/Bag, 50Bags/ Case
RMTBLBS1:1	Blue Light Body Short 1:1	50 pcs/Bag, 50Bags/ Case
RMTRED	Red Dynamic	50 pcs/Bag, 50Bags/ Case
RMTYD	Yellow Dynamic	50 pcs/Bag, 50Bags/ Case

LOUPE

LED HEADLIGHT & BATTERY



With a dedication to serve you best we strive to increase comfort while wearing ORO loupes. We offer loupes with the best resolution across the widest and deepest viewing field. A loupe with exceptional visibility and light weight is a must for long procedures. The hinge system is stable and can be adjusted to satisfy your working position & to support different pupil distance. Every loupe has frame for customers who wear optical prescription glasses.

Rechargeable lithium battery, light weight with long lasting capacity. Easy to carry as well.

SPECIFICATIONS

Magnification	2.5x and 3.5x
Colour	Silver, Red, Blue, Black

LOUPE

LEATHER HEADBAND



This leather headband is light weight and features an adjustable height and width so that a comfortable fit is always assured. The sliding mount allows the user to optimize the position of the optics. Soft fabric cushions add comfort and can be removed for washing, or replacement.

Key Feature

- Best fit with convenience and sturdy construction
- Good light and spot converging
- High intensity and power LED
- Suitable to wear over glasses

HANDPIECE

ORO offers a wide selection of dental hand pieces from high and low speed for a wide variety of dental uses. Our hand pieces are efficient, practical, economical and consistent through selection of very high-quality materials. This provides better sterility and extremely high engineering tolerances lead to cost-efficiency paired with high-quality technology.

We strive to provide the best possible quality service for your small equipment to help you maintain the life of your precision instruments.

GOLD SERIES

- Smooth during use and works for a long time
- Convenient to take off the bur and time saving.
- Internal water spray hole is designed in a way to ensure accurate spray position



SPECIFICATIONS

Code	Product	Rotation Speed Range
GD-S H2	Gold Standard	360-400krpm
GD-M H2	Gold Standard Push Button	360-400krpm
GD-T H2	Gold Super Torque Wrench Head Type	320-360krpm
GD-TM H2	Gold Super Torque Push Button	320-360krpm
GD-LOW01 C	Gold Contra Latch Type	5000-20,000 rpm
GD-LOW01 S	Gold Straight	5000-20,000 rpm

PLATINUM SERIES

- Low noise
- Best performance with integrated head.
- Fantastic optimal access and great accuracy for all treatment scenarios.
- Enjoy the efficiency with comfortable working.



SPECIFICATIONS

Code	Product	Rotation Speed Range
PL-S H2	Platinum Standard Wrench Head Type	380000-420000
PL-M H2	Platinum Standard Push Button	380000-420000
PL-T H2	Platinum Super Torque Wrench Head Type	330000-420000
PL-T3 H2	Platinum Super Torque Push Button	330000-420000
PL-45 H2	Platinum 45° Head	380000-420000

LED HANDPIECE

- Advanced German bearing for maximizing life span of the hand piece.
- Built in with dynamic balancing correction achieving best cutting efficiency, with the clear vision of LED.
- Comes in triple water spray and the LED can withstand 135-degree high pressure sterilization.
- Zero-retraction ensures no cross-contamination

SPECIFICATIONS

Code	Product	Rotation Speed Range
PL-E H2	LED Super torque push button	380000-420000



SURE ENDO

GUTTA PERCHA POINTS



Are meticulously crafted endodontic materials, rolled uniformly from premium Gutta Percha to ensure consistency in size and shape. Compliant with ISO 6877 standards, they guarantee brittleness and radiopacity while remaining flexible and pliable for easy insertion into root canals. Ideal for various endodontic techniques, they enhance irrigant agitation, facilitating effective cleansing and disinfection. Trusted by clinicians, Sure Endo Gutta Percha points offer reliability and performance, essential for successful root canal procedures.

- Standard Points
 - Calibration Marketing
 - Bio G/P Points
 - Resinated G/P Points
- (Patent No. 0961463)

SPECIFICATIONS

Product	Type
Gutta Percha point	2%
	4%
	6%
	Protaper
	Non-Standard

SURE ENDO

WOODEN WEDGES




Paper Points Length Marked are specially designed endodontic tools with a precise length indication system. They are crafted to effectively absorb moisture and facilitate medicament introduction while maintaining consistent quality. Available in various sizes, Sure Endo Paper Points Length Marked provides clinicians with reliable instruments for precise and efficient endodontic treatment.

SPECIFICATIONS


Product	Type
Absorbent paper point	2%
	4%
	6%
	Protaper
	Non-Standard

GUTTA PERCHA & PAPER POINT

ISO STANDARDIZED 0.02




GREATER TAPER 0.04



- 6% / 0.06
- 4% / 0.04
- ISO 2% / 0.02

GREATER TAPER 0.06




Slide Box




4 Comp. Box

PROTAPER SIZE




ACCESSORY SIZE (REGULAR TYPE)



SURE CORD

KNITTED RETRACTION CORD

- Made of 100% cotton, more absorbent than Twist and Braided type.
- Knitted interlocking chains carry significantly greater quantities of hemostatic solution and provide an excellent medium for gel absorption and placement.
- While preparing packing, knitted cord does not entangle in diamond bur.
- Effective for gingival displacement and offers easy packing into sulcus.



CORD CUTTER

- Easy and convenient to cut the cord
- Stainless steel blade inside the cap,
 - Cuts the cord automatically while snapping down the cap
 - Holds the cord inside the cap so the cord does not fall back into the bottle after cutting
 - Protects from cross-contamination
- Patent No. 20-0447858
- Color coded 6 sizes provide quick and easy identification in the sulcus
- Each bottle contains 100 inch (254 cm) of cord
- 000 (XX-Fine), 00 (X-Fine), 0 (Fine), 1 (Medium), 2 (Coarse), 3 (X-Coarse)



SURE CORD PLUS

RETRACTION & HEMOSTASIS

Optimal Cord for Accurate Impression

Impregnated with AICI

- Gentle astringent effect and hemostasis
- Safe - No cardiovascular problem (Epinephrine Free)
 - No blackening of tissue

Smart-Cut




- One touch cutting cap
- Ceramic Cutter
- MM Mark (3cm, 4cm, 5cm) on the bottle

Knitted with Split Microfiber

- Rapid retraction within 1-2 minutes
- Outstanding absorbency
- Stable structure: Easy placement and packing
- No fuzz (lint-free) on the surface of cord

Patent No. 10-2011-0057230, 10-2011-0057229, 20-0447858, 30-0576285

6 Sizes color coded





REACH GLOBAL INDIA PVT. LTD.

Head Office:

Office No. 1 & 2, Survey No. 64/1, Wing-A, First Floor, Sai Raj Regency,
Pimple Saudagar, Pune, State - Maharashtra (India), Pin Code - 411027

Manufacturing Unit:

Plot No 1-22, Gat No. 1064, Chordia Industrial Park, Naigaon, Shirwal,
Taluka-Khandala, Dist. Satara, State-Maharashtra (India) Pin Code - 412801

Phone No: +91 20 2740 0987 | Customer Care No: +91 93707 88850
info@orohealthcare.com | www.orohealthcare.com